	TECHNICAL SUPPORT BULLETIN	Tech. Support Bulletin No. 105
	<u>System Noise Test</u>	Page 1 of 3 Products: Dim Icon
Revision	Description	Date
B	Initial Release – Production (Edited by Kamaljit Singh & Pat McPhail)	12/09/2011

**Subject: System Noise**


**System(s) Affected: Dimension Icon**

**Action Category:**

**Summary: System Noise Test**


1. Verify that the air table is floating properly.
2. Start “tapping in air” experiment.
3. Place a clean bare Silicon wafer on the chuck & place the **Chuck Vacuum switch** to the **ON** position.
4. If using **HOPG** sample on small sample holder, place the chuck Vacuum Switch to the **OFF** position.
5. Place a new **TESPA or OTESPA** probe into a DAFM tip holder.
6. Install the holder onto the head.
7. Center the laser & photo detector to achieve max sum.
8. Tune the probe.
9. Set the following scan parms:

Scan	Channel 1	Channel 2	Feedback
Scan size = 1um	Data Type = Height Sensor	Data Type = Height	Ig = 0.2 Pg = 0.2
Aspect ratio = 1	Data scale = 1nm	Data scale = 1nm	

	TECHNICAL SUPPORT BULLETIN	Tech. Support Bulletin No. 105
	<u>System Noise Test</u>	Page 2 of 3
		Products: Dim Icon
Revision	Description	Date
B	Initial Release – Production (Edited by Kamaljit Singh & Pat McPhail)	12/09/2011

Scan angle = 0°	Line direction = Retrace	Line direction = Retrace	
Scan rate= 1.5hz	Realtime planefit= Line	Realtime planefit= Line	
Samples/Lines = 512	Offline planefit= Full	Offline planefit= Full	
Closed loop = ON			
Z range = 1 um			

10. Engage on the wafer & optimize scan parms to resolve silicon grains.
11. Capture an image named **OL Silicon**.
12. **DO NOT WITHDRAW.**
13. Set the following scan parms:
  - Scan size = 0.1nm**
  - Scan rate = 2.44hz**
  - Sample/Lines = 256/256**
14. Capture an image named **SYSTEM Noise on Silicion at 2.44hz**.
15. **Withdraw the head.**
16. Open the captured image & click on the **Flatten** icon in the top tool bar.
17. Conduct a **1<sup>st</sup> order Flatten** on the entire image.
18. Click on the Roughness icon & verify the **Image Rq result**
19. **ImageRq** should be less than or equal to **.03nm** to pass noise test. If it fails, do a section analysis.

	TECHNICAL SUPPORT BULLETIN	Tech. Support Bulletin No. 105
	<u>System Noise Test</u>	Page 3 of 3
		Products: Dim Icon
<b>Revision</b>	<b>Description</b>	<b>Date</b>
B	Initial Release – Production (Edited by Kamaljit Singh & Pat McPhail)	12/09/2011

20. Conduct a section analysis and draw a horizontal line.
21. Find the peak frequency by pulling the cursor from left side in the histogram.
22. If you see the peak frequency range from 15-40 Hz, it is most likely because of vibration issues. Check you air table and make sure it is floating properly. (40-60 psi).
23. Higher frequencies (100-300Hz) correspond to scanner noise issues. Consult Bruker Tech Support.