



**SHORT ARC LAMP  
OPERATING PROCEDURES**

## SHORT ARC LAMP

Please read these instructions completely before  
operating this equipment.

If there are any questions or problems regarding  
the use of this equipment, please contact:

**ORIEL CORPORATION**  
250 Long Beach Boulevard  
Stratford, CT 06497-0872  
Phone: (203) 377-8282  
Fax: (203) 378-2457

- or -

**ORIEL S.A.R.L.**  
9 Avenue De Laponie  
Z.A. De Courtaboeuf  
91951 Les Ulis Cedex  
France  
Phone: 01-69-07-20-20  
Fax: 01-69-07-23-57

**ORIEL SCIENTIFIC, LTD.**  
1 Mole Business Park  
P.O. Box 31  
Leatherhead  
Surrey KT22 7AU  
England  
Phone: 0372-378822  
Fax: 0372-375-353

- or -

The representative from whom this equipment was purchased.

---

**OPERATING INFORMATION**

---

DO NOT HANDLE OR ATTEMPT TO OPERATE ANY ORIEL SHORT ARC LAMP UNTIL YOU ARE AWARE OF THE HAZARDS INVOLVED.

**I.1 Explosion Hazard**

The pressure inside all Oriell xenon and mercury(xenon) lamps is above atmospheric, even when the lamps are cold. When operating, the internal lamp pressure is several tens of atmospheres. Therefore:

- 1 Extreme care is necessary when handling and mounting a lamp.
- 2 The lamp should only be operated in a safe enclosure, such as an Oriell Lamp Housing.
- 3 The lamp should be mounted in a stress-free manner.

**I.2 Handling a cold lamp:**

For safety, use goggles and heavy gloves when handling a cold lamp. Avoid straining, or twisting the lamp in any way. Never touch a bare lamp with your hands for safety, and to prevent finger oils from being deposited on the lamp. Finger oils will etch the quartz envelope, particularly with UV irradiation. The weakened surface can rupture under the high stress levels in an operation. Prior to operation of the lamp, clean any oil or dirt off with alcohol. Do not touch the quartz envelope after cleaning.

**I.3 Radiation (Light) Hazard**

These arc lamps emit short wavelength ultraviolet radiation and intense visible and infrared radiation. Absorbed radiation (usually ultraviolet and infrared) can cause eye and skin burns and the hazards of ultraviolet light include corneal burns, skin burns and possible mutagenic effects. Therefore:

- 1 Operate the lamp in a manner which confines the radiation to the desired output beam.
- 2 If possible, filter the ultraviolet from the beam to reduce the hazard. If your application prevents ultraviolet beam confinement, then post appropriate warning signs and equip all exposed personnel with appropriate safety wear.

## SHORT ARC LAMP

---

- 3 Never look into the output beam, at any reflection of the beam or at any surface irradiated by the beam without appropriate safety goggles. A "black" surface may be reflective in the ultraviolet or infrared.
- 4 Do not allow an intense, particularly focused, beam to fall on bare skin.
- 5 Do not allow any ultraviolet radiation to fall on unprotected eyes or skin.

### I.4 Electrical Hazard

Transient high voltages are required to start arc lamps, typically 30 kV, and the power supplies provide high currents at moderate voltages. Therefore:

- 1 Ensure that the electrical connections to the lamp are safely enclosed and adequate for the voltages and currents to be used.
- 2 Exercise good electrical safety practice when working with any lamp-power supply combination. If possible, use interlocks to power down the supply if the lamp enclosure is opened.
- 3 Before changing lamps, make sure the power supply is fully discharged. Many arc lamp power supplies have filter or boost capacitors in the output stage; allow time for these to fully discharge before accessing the lamp.

### I.5 Ozone Hazard

Ultraviolet radiation produces ozone, a toxic gas with different physiological impact on individuals and various tolerance levels. Ozone is beneficial in the upper atmosphere but is a ground level pollutant. The American Industrial Hygiene Association publish the following industrial standards for maximum allowable concentrations:

- 0.1 ppm by volume for an 8 hour exposure
- 2 ppm for a 2 hour exposure

Local codes may be stricter. As a guide, 0.015 ppm is barely detectable by smell, and 1 ppm provides a very obnoxious odor. Ozone test kits are available from Oriel. Whatever the code, we do not recommend exposure to ozone. Therefore:

- 1 Use an ozone free lamp unless you require the very short wave ultraviolet.
- 2 Ensure adequate ozone removal by one or more of the following:
  - a Operate in a well ventilated area (large volume).
  - b Use an Oriel Ozone Eater™ for catalytic destruction of ozone from the lamp housing.
  - c Vent the exhaust from the lamp housing, preferably to an ozone destruction system. If you do this, ensure that the lamp cooling is not adversely affected.

### I.6 Lamp Cooling

Arc lamps operate at high internal pressures and temperatures. The envelope and terminals must operate within the design specifications. Terminal temperatures should not exceed 200 °C. The bulb should be uniformly cooled within a range that depends on the lamp type. Proper cooling will ensure the stress on the lamp envelope is that for which the bulb shape was designed. Additionally the operating parameters of mercury arc lamps are sensitive to envelope temperature. Over cooling of these lamps results in high current operation with low output and greatly shortened lamp life. Oriel Lamp Housings are designed for proper operation of these lamps.

### I.7 Lamp Operating Orientation and Electrical Connections

All arc lamps should be operated vertically. For xenon and mercury(xenon) lamps the anode and therefore anode terminal (+ve) must be on top. For mercury lamps the anode and therefore anode terminal should be at the bottom. The +ve lead from the power supply must be connected to the anode.

II

---

**WARRANTY**

---

Oriel's Arc Lamp Warranty is on a pro-rata basis and covers lamps which are operated according to Oriel's instructions. Lamps which are operated incorrectly are not covered. The most common cause of lamp failure is incorrect polarity. Improper cooling, mounting which induces strain and operation of lamps with dirty bulbs are other major causes of premature failure.

We will replace the lamp if it fails within the first 10% of rated life (the average life quoted in the Oriel catalog). If a lamp fails after this initial 10% we will credit you for the unexpired lamp life on a proportional basis. E.g. If the lamp fails at 60% of rated life, we will credit you with 40% of purchase price. Please note that end of rated life is defined as total accumulated operating hours to 50% of initial luminous output. Fill out the attached form and return the lamp to Oriel in accordance with the Return Policy. No claims of any kind will be accepted or processed unless this form is filled out completely and returned with the merchandise in question or under separate cover.

III

---

**RETURNS**

---

Before returning equipment to Oriel for repair/replacement, please call the Customer Service Department. The phone number is (203) 377-8282. You will be given a Return Material Authorization (RMA) number. To save time, have the Purchase Order number used to purchase the equipment available before calling Oriel. Having this number will ensure that your equipment will be properly processed and greatly shorten the time required for the repair/replacement. Equipment returned without a Return Material Authorization (RMA) number may be rejected by the Oriel Receiving Department. Equipment returned under warranty will be returned with no charge for the repair/replacement or shipping.

# SHORT ARC LAMP

Lamp Model No. \_\_\_\_\_ Serial No. \_\_\_\_\_  
Purchased From \_\_\_\_\_ Date \_\_\_\_\_

### Supporting Equipment:

- a) Lamphouse Model \_\_\_\_\_ Mfg. by \_\_\_\_\_
- b) Ignitor Model No. \_\_\_\_\_ Mfg. by \_\_\_\_\_
- c) DC Power Supply Type \_\_\_\_\_  
Model No. \_\_\_\_\_ Mfg. by \_\_\_\_\_

### Operating Conditions:

Accumulated Running Hours \_\_\_\_\_  
Duty Cycle \_\_\_\_\_ No. of Ignitions \_\_\_\_\_  
Voltage \_\_\_\_\_ Current \_\_\_\_\_ Wattage \_\_\_\_\_  
Lamp Axis (Horizontal or Vertical) \_\_\_\_\_  
Cooling (Forced Air, Electrode Cooled) \_\_\_\_\_  
Was the lamp operated in reverse polarity  Y  N  
Conditions causing reject or return \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Additional information \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### For enhanced (boost) current applications only:

Boost pulse duration \_\_\_\_\_  
Boost pulse frequency \_\_\_\_\_ % Boost\* \_\_\_\_\_

$$*\% \text{ Boost} = \frac{i_{\text{max}} - i_{\text{normal}}}{i_{\text{normal}}} \times 100$$

$i_{\text{max}}$  and  $i_{\text{normal}}$  are instantaneous, not r.m.s. values.

Form completed by _____ Title _____
Company _____
Address _____
Phone No. _____ Fax No. _____ Date _____