

# ESCA 8.5 Release Notes

N. Al-Shamma  
8 Apr 1990

This document describes the changes in ESCA between versions 8.4 and 8.5. See Joe Hughes' *ESCA 8.4 Notes* for differences since ESCA 8.31.

## 1. Sensitivity Factors

The "UPDATING SENSITIVITY FACTORS" message and delay are no longer encountered when analyzing data from different X-ray sources.

Sensitivity factors are adjusted using the sensitivity factor exponent to correct for instrumental and escape depth effects. The adjustment is a function of reference energy, which in turn is a function of X-ray source.

Originally, ESCA assumed that the source was always due to the Al K-alpha line and the entire table was adjusted at one time. Then, in a recent version of ESCA, support was added for different sources. This feature was poorly implemented, however, making it very difficult to switch back and forth between data collected with different sources. Every time data from a different source was analyzed, the Scofield Table was read from disk, and all 1300 entries were adjusted for the new reference energy. This would take a minimum of about 10 seconds, and could take minutes on an older computer with a full hard disk.

The software now applies the adjustment only as needed when sensitivity factors are used in the Composition Table or in profile or map analysis. A fully adjusted copy of the factors is used for Survey Analysis. When Survey Analysis is first used after starting the program, or when data from a different source is used, this copy of sensitivity factors is initialized using the the Scofield Table already in memory. No disk access is required.

## 2. Signature Correction

Previously, an extremely non-uniform signature could result in a loss of counts when the signature was applied. To correct this, a "conservation of counts" algorithm was added to the signature correction. This correction is applied at the end of a region or single scan, not during it. The current correction is still applied to displayed data during accumulation. Thus, at the end of a region, the counts may appear to dramatically increase.

The "conservation of counts" algorithm looks at the number of counts in the corrected and uncorrected spectra. It then takes the ratio of counts and applies it to each channel of data so that the total counts in the resulting spectrum is equal to the number of counts measured (i.e., the number of counts in the uncorrected spectrum).

### 3. High Count Rates/Maximum Counts

Accumulation is halted when the total uncorrected counts exceeds 900 million ( $0.9 \times 10^9$ ). Previously, the software could crash if the counts exceeded one billion ( $1.0 \times 10^9$ ), or, the displayed data would be scaled down to one million ( $1.0 \times 10^6$ ), resulting in a confusing display.

The software now checks at the end of every scan, or once a minute for unscanned data, whether the counts exceed 900 million. If they do, the region is considered done, a message is sent to the printer, and accumulation continues to the next region.

### 4. Satellite Removal

An additional satellite removal algorithm has been added which only applies the correction to photoelectron peaks. This is to prevent the application of the removal algorithm to Auger peaks. It is accessed by holding down the shift key when pressing [\*REMU SATELLIT] in the DATA REDUCTION menu. This softkey only appears for non-mono data.

A side-effect of this algorithm is that the Composition Table and Survey Analysis tables are cleared when it is used. This is not considered serious since the satellite removal will normally be done before any other data analysis.

The new algorithm uses the Find Peaks routine (in the Survey Analysis menu), then applies the satellite removal algorithm to the channels within the baseline of photoelectron peaks. An additional test is added for loss peaks, since they can be put in the same category as Auger peaks by the Find Peaks routine. To be considered a loss peak, a peak must have photoelectron peak of at least twice its height which is -10 to -28 eV away.

### 5. Aperture

The new aperture is supported. This aperture has two leaves. Aperture setting one moves the inner leaf into position. Aperture setting two moves both leaves into position. The software applies a one second delay whenever one of the leaves is moved. Each leaf is moved individually.

## 6. Ultraviolet Photoelectron Spectroscopy

A new option was added to the Spot setting in single scan setup and the element table. It is accessed by typing "UPS" for the Spot. When entered, the X-ray source will be turned off during accumulation, and UPS will be printed in printouts and plots. "UPS" may also be entered in lower case, or as just the first one or two characters. For example, "u" or "UP" will also be interpreted as the "UPS" source.

## 7. DIAGNOSE


Two changes were made:

- a. The dispersion constant for Res 1 is just set to half of that for Res 2. The partially developed algorithm which was erroneously shipped with the product has been removed.
- b. The instrument constants for the V1 curves are now determined by DIAGNOSE. A new softkey, [Calc V1 curves], appears on the top menu if the 8701B power supply is detected.

## 8. Composition Table Printout

The file description and collection date are now printed out on the long Composition Table printout. (No change was made to the short printout.)

		REVISIONS			
LTR	DESCRIPTION	EO	DATE	APPROVED	
B	TYPICAL CHGD SEE ECD	392	6-9-84	J.F. W...	
C	MISC DESIGNATOR CHGS	B23	5/12/86	J. FOX	



RESOLUTION	NOMINAL V OUT
1	25 5
2	50 25
3	100 50
4	150 100

6-8-90  
RHunt


060050-UPS/

NOTE:  
 1.  $\nabla$  F DENOTES FLOATING COMMON  
 2. PC BOARD MARKED BY UPPER RIGHT CARD EXTRACTOR IS GREEN

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION WHICH SHALL NOT BE REPRODUCED, TRANSFERRED TO OTHER DOCUMENTS, DISCLOSED TO OTHERS, USED TO MANUFACTURE, OR ANY OTHER PURPOSE WITHOUT PRIOR WRITTEN PERMISSION OF SURFACE SCIENCE LABORATORIES.

QTY		ITEM		DESCRIPTION	
CONTRACT NO.					
APPROVALS	DATE	SURFACE SCIENCE LABORATORIES			
DRAWN GDI	1-15-84	TITLE SCHEMATIC, PCB			
CHECKED PEL	6/17/84	SPHERE SUPPLY, SPECTROMETER P.S.			
MATERIAL		SIZE DRAWING NUMBER			
FINISH		D SCH-060050-00 C			
NEXT ASSEMBLY		USED ON		SCALE NONE	
APPLICATION		DO NOT SCALE DRAWING		SHEET 1 OF 1	

REVISIONS					
LTR	DESCRIPTION	EO	DATE	APPROVED	
C	TRACES CHANGED - C11-118 W/ 100 FT		291/431	6-9-84	RLI/4/44
D	ADDED 8 .01 CAPS C17-C24		510	10-5-84	MTE/1/1/1/1
E	ADDED NOTES 3 & 4 UPDATE BOM RELEASED BACK TO PRODUCTION		510	10-1-86	HG



RESOLUTION	NOMINAL V OUT	BE
4	500 334	100
3	334 167	50
2	167 83	25
1	83 16	5


060060-UPS

- NOTES: UNLESS OTHERWISE SPECIFIED
- ↓ F DENOTES FLOATING COMMON
  - PC BOARD MARKED BY UPPER RIGHT CARD EXTRACTOR IS BLUE
  - CAPACITANCE VALUES ARE IN MICROFARADS
  - RESISTORS ARE .125 WATT ±1%

QTY		ITEM		DESCRIPTION	
CONTRACT NO.					
APPROVALS		DATE		TITLE	
DRAWN		7-21-85		SCHEMATIC, PCB, OFFSET SUPPLY, SPECTROMETER POWER SUPPLY	
CHECKED				SIZE	
FINISH				DRAWING NUMBER	
I-8701				SCH-060060-00	
NEXT ASSEMBLY		USED ON		REV	
APPLICATION		DO NOT SCALE DRAWING		E	
		8701		SCALE	
				SHEET 1 OF 1	

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION WHICH SHALL NOT BE REPRODUCED, TRANSFERRED TO OTHER DOCUMENTS, DISCLOSED TO OTHERS, USED TO MANUFACTURE OR ANY OTHER PURPOSE WITHOUT PRIOR WRITTEN PERMISSION OF SURFACE SCIENCE LABORATORIES.

A	R8, 9, 10, R11 WAS 100K	3/2/83	5-30-81	WLB/1
B	ADDED .01 CAPS, C13 THRU C20	5/11	10-5-84	MTE 1/1/85
C	ADDED NOTES 1 & 2 AND VOLTAGE TO CAPS RELEASED BACK TO PRODUCTION.	5/11	9-30-86	H. G. 7/1/86



NOTES: UNLESS OTHERWISE SPECIFIED  
 1. CAPACITANCE VALUES ARE IN MICROFARADS.  
 2. RESISTORS ARE .125 W ± 1%


### 060080-UPS

RESOLUTION	OUTPUT I
4	136 mA ± 50
3	110 ± 50
2	90 ± 50
1	70 ± 50

6-8-90  
R Hunt

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION WHICH SHALL NOT BE REPRODUCED, TRANSFERRED TO OTHER DOCUMENTS, DISCLOSED TO OTHERS, USED TO MANUFACTURE, OR ANY OTHER PURPOSE WITHOUT PRIOR WRITTEN PERMISSION OF SURFACE SCIENCE / KEYEX CORP.

GDI	2/22/83	SCHMATIC, PCB, MAGNET SUPPLY, SPECTROMETER P.S.	SIZE D	REV C
		870		




REVISIONS				
REV	DESCRIPTION	ED	DATE	APPROVED
C	TRACES... CHANGED... 291/431	67-84	RL	
D	ADDED B... OF CAPS... 510	10-5-85	MTB	
E	ADDED... 510	10-1-86	MC	

RESOLUTION	NOMINAL V OUT	BE
4	500 334	100
3	334 167	50
2	167 83	25
1	83 16	5


060060-UPS

- NOTES: UNLESS OTHERWISE SPECIFIED
1.  $\nabla$  DENOTES FLOATING COMMON
  2. PC BOARD MARKED BY UPPER RIGHT CARD EXTRACTOR IS BLUE
  3. CAPACITANCE VALUES ARE IN MICROFARADS
  4. RESISTORS ARE .125 WATT  $\pm$  1%

DATE: 7-27-01		DRAWN: JLM		DATE: 7-27-01	
CHECKED: [ ]		APPROVED: [ ]		DATE: [ ]	
MATERIAL: [ ]		TITLE: SCHEMATIC, PCB, OFFSET SUPPLY, DELETED NET LBL, POWER SUPPLY		SCALE: [ ]	
PARTS LIST: [ ]		DRAWING NUMBER: D SCH-060060-00		REV: E	
APPLICATION: [ ]		JOB CODE: 8701		SHEET 1 OF 1	



- A RB, 10, 11 WAS 100 K
- B ADDED  $\Delta$  .01 CAPS, C13 THRU C20
- C ADDED INPUT 1 & 2 AND VOLTAGE TO CAPS. THIS WOULD HAVE TO BE REWORKED



NOTES: UNLESS OTHERWISE SPECIFIED  
 1. CAPACITANCE VALUES ARE IN MICROFARADS.  
 2. RESISTORS ARE .125 W = 1%

**060080-UPS**


RESOLUTION	OUTPUT I
4	136 mA ± 50
3	110 ± 50
2	90 ± 50
1	70 ± 50

6-8-90  
 RHUNT

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION WHICH SHALL NOT BE REPRODUCED, TRANSMITTED TO OTHERS, DISCLOSED TO OTHERS, USED TO MANUFACTURE OR ANY OTHER PURPOSES WITHOUT PRIOR WRITTEN PERMISSION OF SURFACE SCIENCE LABORATORIES.

GO1	2/22/83	SCHEMATIC, PCB, MAGNET SUPPLY, SPECTROMETER, 25, 2500	SIZE B	SURFACE SCIENCE LABORATORIES
		8701	SIZE B	SCH-060080-001C





LTN	DESCRIPTION	REV	DATE	APPROVED
B1	TYPICAL CHGD SERA CSD	392	6-7-80	PELLUM
C	MISC DESIGNATOR CHGS	B23	5/16/86	OFFER

RESOLUTION	NOMINAL V OUT
1	25 5
2	50 25
3	100 50
4	150 100

6-8-90  
R Hunt  
060050-UPS/

- NOTE:
1.  $\nabla$  DENOTES FLOATING COMMON
  2. PC BOARD MARKED BY UPPER RIGHT CARD EXTRACTOR IS GREEN

QTY	ITEM	DESCRIPTION

CONTRACT NO	DATE	APPROVALS
DESIGNED BY	DATE	APPROVALS
CHECKED BY	DATE	APPROVALS
DATE	SCALE	SCALE

SURFACE SCIENCE LABORATORIES  
 TITLE: SCHEMATIC, PCB  
 SPHERE SUPPLY,  
 5 SPECIMEN  
 SIZE: D 5CH060050-00 C  
 DRAWING NUMBER: 8701  
 SCALE: NONE  
 SHEET: 1 OF 1