


Degassing the X-Ray filaments

1. Turn off X-Ray Controller (9600), the Glassman High Voltage Power supply.
2. Check
 - a. that the Hawk water circulator, is on and ready.
 - b. that the cryo is on and down to temperature and that the gates to the cryo are open.
 - c. that the vacuum is below 5×10^{-8} Torr.
 - d. that gate 5 is closed
3. Turn on Glassman high voltage and the 9603 X-ray gun spot size controller
4. Make sure the "Interlock OK" LED comes on otherwise you cannot get current from gun.
5. Turn "**ramp**" knob on spot size controller to slowest, **fully clock wise** set tge "**service**" switch **down**
6. Press "**start filament**" on 9600 followed by "**HVon**" on Glassman It will go through 2 cycles of ramping voltages?
7. Ramp 1:
 - a. The front panel "Fil on" LED should turn on and the Panel meter should read I FIL mode and go from 0.6 to ~ 1.2 A. Inside the 9600 XRay controller the "Fil" (05) LED should be on and LED 06 will come on
8. Ramp 2, 2KV:
 - a. Meter V2KV go to 2.3 KV
 - b. LEDs 05 off, 06 off, 04 (2KV) on, 07 will come on
9. Watch that the Pressure does not go up too much.
10. Set the "**service**" switch **UP**, keeping "**ramp**" knob on slowest setting (**c.w.**) and changing "**stand by**" to "**operate**". The Glassman should now read about 2 KV and the Voltage control should be active. This will take ~ 8 h, watch pressure.
11. Once 10 kV is reached, degas anode by starting X-ray gun on **100 micron** spot and slowly increasing the spot size until largest spot does not raise pressure above 2×10^{-8} Torr
12. Turn off X-rays
13. Turn on spectrometer boxes (top 3, spectrometer power supply, flood gun, memory interface)



14. Open software, check "X-ray gun operate" in ESCA control panel
15. Switch spot size controller from "**manual**" to "**computer**"
16. Degas flood gun - check "**Flood Gun**" box in ESCA control panel and increase energy to 5 eV, watch pressure and wait for at least 1 hour
17. Turn down flood gun energy, uncheck **flood gun** box, and close ESCA control panel
18. Instrument is ready to use