

## Restart of MProbe For GLA's only


### Mprobe Restart After Vacuum Loss

- 1) Turn off all electronics, (2505 Memory Interface, Spectrometer Power Supply, Glassman High voltage, Service Physics X-Ray Gun Controller, nti sputter gun controller)
- 2) Turn off alarm on interlock box near bottom of rack and set all three bypass switches on box **up**.
- 3) Push **reset** button on interlock box
- 4) Press **reset** on back of control box over XPS
- 5) If the Neslab circulator is off **Restart** it.
- 6) Check that the transfer arm is full withdrawn.
- 7) **If the Cryo's are up to a temperature >80 K or the pressure is higher then  $1 \times 10^{-7}$  Torr.**
  - a) Turn **on** interlock override on control box.
  - b) Turn **off** cryo-pump compressor
  - c) If the Turbo pump is on
    - i) Shut gate 2 and turn Turbo **off**.
    - ii) Bring chambers up to atmosphere, Open gate 1, 5, and 3, b
  - d) Wait for the cryo's to warm up  $T(\text{cryo}) > 270\text{K}$  (will take a few hours);
  - e) Open gates to cryo 0 and 6.
  - f) Shut Nitrogen, gate 3.
  - g) Pump the system, Open gate 2, Turn **on** the turbo pump.
  - h) Wait till the vacuum is  $< 1 \times 10^{-6}$  T, this will take 1 to 2 hours or more.
- 8) Turn **on** the cryos (you may need to turn the cryo switch off and then on).
- 9) Make sure that the cryo pump is working
- 10) After the cryo temperature is down below 25K (about 1 to 3 hours) **close** gate valves 1 and 5 (gates 0, 2, 5, and 6 are open).
- 11) On the interlock box set the alarm and three bypass switches to **down** (on)
- 12) Turn **off** interlock override.
- 13) Continue to next section.

### Restarting M-Probe X-Ray Gun

- 1) Check that the vacuum system is OK
- 2) Gates 0, 2, and 6 should be open with all others closed
- 3) Turn **off** interlock override on back of control box over XPS.
- 4) Make sure the system power switch on the back is **up**
- 5) Turn **on** X-ray gun boxes
  - a) Turn on Glassman high voltage and press HV on,
  - b) Turn on 9603 X-ray gun spot size controller

- c) Make sure the “**Interlock OK**” LED comes on otherwise you cannot get current from gun.
- 6) If the vacuum has gotten above  $>1 \times 10^{-7}$  Torr or system was vented:
- Turn "ramp" knob on spot size controller to **slowest**, fully clock wise with the service switch **up**. Note that the switch in the up position runs the ramps below about 100 times **slower** than when it is in the **down** position.
  - Press "**start filament**" on 9600 followed by "**HV on**" on Glassman It will go through 2 cycles of ramping voltages? If the pressure does not rise after 2 minutes you can flip the switch **down** and watch the pressure. If it increases flip the switch **up** and wait till the pressure goes down. You can repeat this till the ramps are complete,
    - Ramp 1, I FIL:
      - The front panel “**Fil on**” LED should turn on and the Panel meter should read I FIL mode and goes from 0.6 to  $\sim 1.2$  A.
      - Inside the 9600 X-Ray controller the “**Fil**” (05) LED should be on and LED 06 will come on
    - Ramp 2, V 2KV:
      - Meter V 2KV goes to 2.3 KV
      - LEDs 05 off, 06 off, 04 (2KV) on, 07 will come on
  - Watch that the Pressure does not go up too much if it does flip the switch **up**.
  - Set the "service" switch **UP**, keeping "ramp" knob on **slowest** setting (c.w.) and changing "stand by" to "**operate**" (Make sure the “**Interlock OK**” LED comes on otherwise you cannot get current from gun. It takes  $\sim 8$  h, watch pressure) The Glassman should show 1 or 2 KV of high voltage but almost no current.
- 7) If the vacuum did not get above  $<1 \times 10^{-7}$  Torr:
- Turn "ramp" knob on spot size controller to **fastest**, counter clock wise (service switch down)
  - Press "**start filament**" followed by "**HV on**" and for the filament and 2KV ramp to complet (see 6b above) This should take a few minutes.
  - Watch the pressure if it increases significantly go to 6 above
  - Set the "service" switch **up** and turn "ramp" knob to **fastest**, then press "**operate**" (takes 2 h, watch pressure)
- 8) Once 10kV reached, degas anode by starting X-ray gun on 100 micron spot and slowly **increasing the spot size until largest spot does not raise pressure above  $2 \times 10^{-8}$  Torr**
- 9) Turn off X-rays
- 10) Turn on spectrometer boxes (top 3, spectrometer power supply, flood gun, memory interface)



- 11) Open software, check "**X-ray gun operate**" in ESCA control panel
- 12) Switch spot size controller from "manual" to "**computer**"
- 13) Degas flood gun - check "**Flood Gun**" box in ESCA control panel and increase energy to 5 eV, watch pressure and wait for at least 1 hour
- 14) Turn down flood gun energy, **uncheck flood gun** box, and close ESCA control panel
- 15) Instrument is ready to use

### **XPS maintenance schedule Kratos and MProbe**

- **Every week:**
  - Refill water on Affinity and Neslab circulators
  - Check water level on Hawk chiller
  - Check N<sub>2</sub> level on HREELS and order new cylinders
  - Replenish gloves, IPA, and KimWipe
- **Every 6 months:**
  - Change M-Probe pump oil (last 08/2019)
  - Change M-Probe anode (3 3/8" Cu gasket, Al/Cu anode, 6 Au 1/8" screws PRT-6520-003) (last 08/2019)
  - Change M-Probe Hawk circulator deionizer (Thermo Sci combined DEM/OXY cartridge D8809) and water filter (Hytrex cartridge filter PRT-6530-003) (last 08/2019)
  - Change the water and air filters for the MProbe located in the overhead rack.(changed 08/2019)
  - Change the oil in the MProbe Fore pump.
  - Change Kratos Affinity circulator filter (Pentek filtration polydepth filter cartridge PD-1-934) (last 2013-09-25)
- **Every 36 months (last 2013-02-17):**
  - Check and change Kratos deionizer cartridge (84-789), or when deionizer stays on
  - Change house water filter
- **Every 60 months:**
  - Change tip seals on Kratos (last 2013-04-19) and EELS (last 2010-07-01) scroll pumps
  - Change house air filter (Motor guard M-723 filter element) (last 2010-07-01)

### **Contacts for M-Probe: Service Physics (Bend, Or) (541) 318-8688**

**Zack Mehl** [Zach@sphysics.com](mailto:Zach@sphysics.com) - New service Engineer, Mike's replacement

**Bob Chaney:** bob@sphysics.com - the boss, good contact email or phone for software or in general x1

**Barbara Siordia:** barb@sphysics.com - sends packages x2

**Ruth Chaney:** ruth@sphysics.com - Bob's wife, does POs 541 322 9405