## **Glove Box Standard Operating Procedures:**

## **Inserting objects into box:**

- 1. Check to se if anyone else is using box.
- 2. Check to see the inner port door is closed.
- 3. Check to see if port is under vacuum.
- 4. Shut off vacuum (VAC valve should be almost perpendicular to pipe and parallel to port
- 5. Turn on  $N_2(g)$  (small green valve, should now be parallel o pipe)
- 6. Turn off  $N_2(g)$  when round gauge above the VAC reads approximately -1 psi (below 0).
- 7. Open outer port door. This should require force if a slight vacuum is maintained in step 4.
- 8. Insert your stuff into port and close the door.
- 9. Pull VAC valve parallel to the pipe to evacuate port.
- 10. Watch pressure gauge go down below -30"  $H_2O$ .
- 11. Backfill with  $N_2(g)$  again (Steps 3 and 4).
- 12. Repeat steps 9 11.
- 13. Repeat steps 9 11.
- 14. Backfill with  $N_2(g)$  until pressure is equalized between box and port (pressure can go to 0 psi).
- 15. Open inner port door and remove your stuff and close port door.
- 16. Put the port under vacuum again.

## Removing objects from the box:

- 1. Check to see that outer door is closed.
- 2. Check to see if port is under vacuum.
  - a. If it is do steps 4, 5 and 14 above.
- 3. If port is not under vacuum check to see if it was open to air last.
  - a. If port contains air  $(O_2)$  then repeat steps 9 through 13 above.
- 4. If port now contains  $N_2(g)$  open inner port door and place your objects inside.
- 5. Close inner port door.
- 6. Open outer port door and remove objects.
- 7. Pull VAC valve parallel to the pipe to evacuate port.