

## Gamry Instruments' .DTA File Type

### Introduction

Our Framework software is the controlling application for all Gamry Instruments potentiostats. When data are stored for later analysis and plotting, they are saved as a .DTA file type. This Technical Note describes more about the .DTA file type, and how to transfer data into other software packages such as spreadsheet programs.

### The .DTA File Type

Gamry Instruments uses a proprietary file type with the extension .DTA for storing Framework™ software data and electrochemical impedance spectroscopy data generated in Echem Analyst™ software. (For a general discussion of file types in our software, see our [“Software”](#) webpage.)

A .DTA file contains only tab-separated raw text, and can be opened in any text editor or spreadsheet program.

### Column Headers in .DTA Files

Gamry Instruments' .DTA files include headers. The following abbreviations are found in .DTA file headers:


Pt	Data point number, starting with point 0 (not point 1)
T	Time since the start of the experiment
Vf	Measured cell voltage
Im	Measured cell current
Vu	Uncompensated voltage
Sig	Voltage from the signal generator entering the current amplifier
Ach	Voltage measured using the A/D input
IERange	Which current range was used for measurement
Over	A code indicating various kinds of error conditions

### Copy and Paste Method in Echem Analyst™

You can copy and paste a section or entire results sheet from Echem Analyst.



The **Copy to Clipboard** button (see right) is on the toolbar of the tab showing the graph when your data file is displayed in Echem Analyst. This button gives you the choice to copy a Bitmap, Metafile, or raw data directly to the Windows® Clipboard. You can then paste the data into any other program. The **Copy to Clipboard** button opens a drop-down menu of choices (Fig. 1).



**Figure 1.** Copy to Clipboard drop-down menu.

- The **As a Bitmap** option copies the data as displayed in Echem Analyst, meaning that the data will be a replica of the current view of the data complete with any overlays, fits, or modifications of the visual style of the plot.
- Using the **As Text (data only)** button copies the x and y values of the currently displayed data, fits, and overlays as a table that you can directly paste into a cell-based spreadsheet application. By overlaying multiple data files you can use the Copy As Text feature to batch-export a number of data files all at once.

You can copy data to the clipboard from any of the **Analysis** tabs (e.g., Fig. 2) that are created when a fit is performed. Click in the cell (or click and drag to highlight multiple cells; see Fig. 3), then right-click and select **Copy**. You can now paste the data into any cell-based spreadsheet program.

Parameter	Value	± Error	Units
R1	24.88e-3	444.4e-6	ohm
R2	23.75e-3	1.106e-3	ohm
Yo4	1.813	305.8e-3	S*s^a
a5	573.3e-3	31.91e-3	
R6	48.11e-3	6.692e-3	ohm
Yo7	246.8	38.65	S*s^a
a8	685.6e-3	48.22e-3	
Yo9	1.856e3	154.1	S*s^a
a10	1.000	11.78e-3	
L10	757.4e-9	12.87e-9	H
Goodness of Fit	480.6e-6		
13p_pstaticEIS_100uHz.DTA			

**Figure 2.** A typical **Analysis** tab with data created after a fit of a model.

Parameter	Value	± Error	Units
R1	24.88e-3	444.4e-6	ohm
R2	23.75e-3	1.106e-3	ohm
Yo4	1.813	305.8e-3	S*s^a
a5	573.3e-3	31.91e-3	
R6	48.11e-3	6.692e-3	ohm
Yo7	246.8	38.65	S*s^a
a8	685.6e-3	48.22e-3	
Yo9	1.856e3	154.1	S*s^a
a10	1.000	11.78e-3	
L10	757.4e-9	12.87e-9	H
Goodness of Fit	480.6e-6		
13p_pstaticEIS_100uHz.DTA			


**Figure 3.** Data in the **Analysis** tab are selected, with a small menu that appears with choices to right-click.

### Directly Opening the File in Your Spreadsheet

Many customers extract their .DTA data via the copy-and-paste method described above. But there is another possibility: directly opening the .DTA file in your spreadsheet program. We sometimes find that opening a .DTA file in Excel® or another

program works more reliably than copying and pasting.

For example, Fig. 4 is a .DTA file for an EIS experiment opened directly in Excel®:



**Figure 4.** EIS file in .DTA format opened directly in the Excel® spreadsheet program.

### Conclusion

The .DTA file type for Gamry Instruments' software does not have to be treated as a mysterious black box. To transfer your data to other programs, you can simply select, copy, and paste; or you can open the file in your spreadsheet application.

Application Note Rev. 1.0 1/2/2019 © Copyright 2018 Gamry Instruments, Inc. Framework and Echem Analyst are trademarks of Gamry Instruments, Inc. Excel® and Windows® are registered trademarks of Microsoft Corporation.

