

**BigCatch – EN Series Eyepiece Software  
Quick Installation Guide Windows XP or 2000**

1. Copy the Software in software folder to CD or Thumb drive
2. Turn on your PC and place the CD into the CD drive. **Do not connect the USB cable of microscope to PC.**
3. Go to the drive and double click Autorun.
4. Click on Scopephoto and install
5. Plug in the EM-500M to the USB cable and cable into PC.
6. A “Found new hardware wizard” will appear.
7. Tell it you want to search a folder for the driver and Browse to “folder 500USB(2.0).
8. Click next.
9. Click “DirectX9.0” in the Scopephoto disk to install the camera’s Video program.
10. Restart your computer.
11. Click on Scopephoto icon on desktop to run program.
12. **Make sure you close the program before disconnecting the microscopy eyepiece camera from the PC.**