

November 11, 1998

Repair of the Neslab Coolflow System II

Videl from McCullum Air inc. (office phone (626) 795-6380; worked from 7:30 AM to 1:45 PM) disconnected the system from the wall plumbing, removed the bad motor and attached pump and replaced it with the motor and pump from the "extra" leaky Neslab System II Alan had obtained from somewhere. The motor did not start initially, but after some coaxing did, but became hot almost immediately. Oil was added and the pump restarted. The motor did not become as warm, but was still very warm. Videl removed the pump and replaced it with one originally in the unit because that one turned more easily. The motor again needed some coaxing to get it to start, and once started became very hot. Tony from the Biology refrigeration repair shop says this can be normal with this type of motor. The Neslab was left running to see if it would overheat. Meanwhile a new motor was ordered through Tony to replace the bad motor.


ROOM COOLING Loop

NESLAB POWER-FAIL SHUTDOWN.

6-95

NESLAB RECIRCULATOR PUTS COLD WATER INTO OVERHEAD PIPES FOR COOLING LAB EQUIPMENT. WHEN THE NESLAB SHUTSDOWN THIS OVERHEAD WATER DRAINS BACK INTO THE NESLAB TANK CAUSING IT TO OVERFLOW.

A CHECK VALVE HAS BEEN PLACED IN THE "UP" PIPE (SUPPLY-SIDE). AN ELECTRIC SOLENOID VALVE HAS BEEN PLACED IN THE "DOWN" RUN (RETURN SIDE). THE CIRCUIT BELOW CLOSES THIS SOLENOID VALVE WHEN THE NESLAB SHUTS OFF.



7-9-96

"POWER OFF LATCH"

9-96

NESLAB MODIFICATION

AS ORIGINALLY WIRED, THE NESLAB WOULD LATCH ~~ON~~ ^{OFF} FOLLOWING A POWER FAILURE. IT HAS BEEN MODIFIED TO RESTART WHEN POWER RETURNS. THIS WAS DONE AT THE TERMINALS OF THE FRONT PANEL POWER ON/OFF TOGGLE SWITCH. THE SWITCH HAS TWO INDEPENDENT CIRCUITS. CIRCUIT "ONE" (INBOARD) IS MAIN POWER ON/OFF. CIRCUIT "TWO" (OUTBOARD) WAS THE MOMENTARY CONTACT USED FOR INITIAL ~~ON~~ ACTIVATION OF THE POWER-ON LATCH. THE TWO LEADS ~~W~~ ATTACHED TO SWITCH CIRCUIT "TWO" WERE SIMPLY TIED TOGETHER BY MOVING THE BLACK LEAD UP TO JOIN THE YELLOW AT THE UPPER TERMINAL. CAUTION - THIS PROBABLY DISABLES THE SAFETY CIRCUITS, SUCH AS, LOW WATER.

NESLAB BI 023

IN-ROOM RECIRCULATOR

PUMP SHUTDOWN - CAUSED

1/95 OVERFLOW OF RESERVOIR ONTO FLOOR. REPLACED MOST OF THE WATER IN SYSTEM W/ DI WATER.

3/95 REPLACED EXTERNAL GREEN GARDEN HOSES.

INSTALLED SOLENOID VALVE ON ^{INSTRUMENT} NO WATER RETURN LINE AND ONE-WAY CHECK VALVE ON INSTRUMENT SUPPLY SIDE, INTERLOCK TO NESLAB OPERATION. WHEN NESLAB TURNS OFF THE SOLENOID WILL CLOSE TRAPPING WATER IN OVERHEAD PLUMBING.

WATER IS TRAPPED IN THE OVERHEAD PLUMBING WHEN THE NESLAB TURNS OFF FOR ANY REASON. THERE IS A CHECK VALVE ON THE SUPPLY SIDE (GOING UP) AND AN ELECTRIC SOLENOID VALVE ON THE RETURN SIDE (LINE COMING DOWN). THESE VALVES SHOULD KEEP THE NESLAB RESERVOIR FROM OVERFLOWING DURING A POWER OUTAGE.

7-1-95 QUALITY OF WATER $\approx 77 \mu S$. FLUSHED DI ($C = 0.0 \mu S$) THROUGH SYSTEM. AFTER FLUSH $C \sim 5-10 \mu S$. ADDED 10 mL OF VWR CLEAR BATH ALGECIDE TO REDUCE GROWTH INSIDE CLEAR PLASTIC HOSES. CONDUCTIVITY MEASURED WITH H_2O PURE HAND HELD DIPPER UNIT IN TOOLBOX.

5-8-96 CONDUCTIVITY WAS $40 \mu S$. FLUSHED SYSTEM WITH DI HOUSE WATER TO $C \sim 3 \mu S$. (MISTAKE - OPENED CROSS CONNECTS OVERHEAD BETWEEN SUPPLY & RETURN INTRODUCING LARGE PIECES OF BROWN SCALE INTO SYSTEM. NOT GOOD FOR TURBO PUMP COOLING LINES.) ADDED 20ml ALGECIDE (ALSO CAUSES FOAMING AND STRANGE SURFACE TENSION PROPERTIES ON WETTED SURFACES)

"OMNIFILTER!"

5-16-96 ADDED BABY BLUE PARTICLE FILTER OVERHEAD AT XPS FEEDPOINT. 25 μ CARTRIDGE INSTALLED. CARTRIDGES ARE STANDARD 9 $\frac{3}{4}$ " ELEMENTS (AVAILABLE FROM McMASTER-CARR). 20 μ PAPER FILTER

6-10-96 TOTAL FLOW TO M-PROBE $\sim 2.4 \text{ gal/min}$.

7-10-96 OILED BOTH ELECTRIC PUMP MOTORS.

OMNIFILTER ELEMENT WAS OLIVE-GREEN; REPLACED. BEWARE OF WATERFALL ON OPENING FILTER CANNISTER. 15-20 DROPS OF SAE 10 @ EACH BEARING.

[CLOSING THE WALL MOUNTED RETURN LINE SOLENOID LINE CAN IMPROVE FLUSHING BY NOT RECIRCULATING OLD WATER]

4-9-97 REPLACE OMNIFILTER ELEMENT [OLD ONE WAS QUITE GREEN & ALIVE]. HI CAPACITY PLEATED PAPER (20μM) ELEMENT INSTALLED, WATER SAMPLED FROM WEST SIDE SYSTEM DRAIN WAS 91μS. TOTAL FLOW WAS ~ 2.7 gal/min. ^{FLUSH} SUPER DOSE OF ALGECIDE RUN THROUGH SYSTEM. FLUSH AGAIN w/ DI. WATER TO ~ 5μS.

1-14-98 OMNIFILTER ELEMENT WAS GREEN. WATER IN NESLAB TANK WAS CLEAR 46μS. INSTALLED RESIN-BONDED CELLULOSE FILTER TEMPORARILY AND ADDED 1L OF HOUSEHOLD BLEACH TO SYSTEM FOR 3 HOUR (C >> 100μS). FLUSH @ ~ 2 g/min FOR 1 HOUR. LUBRICATED BOTH PUMP MOTORS. INSTALLED 20μM PLEATED PAPER FILTER. FLUSHED AGAIN NEXT MORNING TO 30μS CONDUCTIVITY. CHLORINE TREATMENT CAN DISLODGE QUITE A BIT OF GREEN JUNK TAKING SEVERAL HOURS TO FILTER OUT.

8-25-98 THE RETURN LINE ANTI-DRAIN VALVE WAS STUCK IN THE OPEN POSITION. REPLACED IT WITH A NEW ONE. SPARES (2) AVAILABLE IN LAB.

11-11-98 THE Neslab Coolflow System # was repaired (see note 3 pages back). After changing the pump some foaming was apparent in the tank. A DI WATER exchange was performed for 1 h with the Return line closed (~2.8 g/min) most of the time (40 min), and with the return line open (~2 g/min) intermittently (20 min). THE CONDUCTIVITY meter needs batteries so did not measure the conductivity.

11-12-98 conductivity ~ 70μS

12-98 FILTER REPLACED IN BLUE CANISTER (By Nick)

liquid or gas by forcing them through a porous medium. The particles are measured or classified in *microns*. The micron equals .000039 of an inch. Microns are some linear equivalents and compare the relative size of particles.

1 inch.....	25.4 millimeters.....	25,400 microns
1 millimeter.....	0.0394 inches.....	1,000 microns
1 centimeter.....	1/25.4 of an inch.....	10,000 microns
1 micron.....	3.94 x 10 ⁻⁴ inch.....	0.000039 inches

Relative Sizes

Diameter of a human hair.....	70 microns
Lower limit of visibility (naked eye).....	40 microns
White blood cells.....	25 microns
Red blood cells.....	8 microns
Bacteria (cocci).....	2 microns

General Rating

100 microns.....	Very coarse
50 microns.....	Coarse
20 microns.....	Medium
10 microns.....	Fine/medium
5 microns.....	Fine
1 micron.....	Very fine
.1 micron.....	Submicron

Multipurpose String Wound Filter Cartridges

Filter cartridges offer a much greater area of filtering material and can hold a larger volume of dirt particles. Carbon filters reduce the taste and odor in drinking water. Membrane and submicron filters are used in high efficiency applications. The cartridges listed fit most 4", 10", 20", and 30" filter housings. The ends are Double Open Ended (D.O.E.) unless otherwise indicated.

Wound Polypropylene

- Multi-purpose for organic and petroleum oils
- Non-FDA grade • Polypropylene core
- Maximum temperature is 130° F • 2 1/2" OD

Full-Body Wound Polypropylene

- Water and general filtering
- Have extra-large filtering capacity
- Fit large diameter housings • Polypropylene core
- Maximum temperature is 130° F • 4 1/2" OD

Submicronic Wound Cotton

- Pre-filter for membranes, high purity solvents, and critical water
- Polypropylene core • Maximum temperature is 120° F • 2 1/2" OD

Wound Acetate	Wound Rayon	Wound Cotton	Wound Polypropylene	Full-Body Wound Polypropylene
4411K13.....\$2.93	4411K14.....\$1.82	4411K31.....\$3.38	4411K71.....\$1.89	4443SK11.....\$4.11
4411K21.....\$3.91	4411K81.....\$2.40	4411K32.....\$3.44	4411K72.....\$1.96	4443SK12.....\$4.18
4411K22.....\$3.96	4411K82.....\$2.47	4411K33.....\$3.49	4411K73.....\$2.03	4443SK13.....\$4.25
4411K23.....\$4.18	4411K83.....\$2.67	4411K34.....\$3.52	4411K74.....\$2.08	4443SK14.....\$4.32
4411K24.....\$4.53	4411K84.....\$2.93	4411K35.....\$3.61	4411K75.....\$2.14	4443SK15.....\$4.39
4411K25.....\$4.80	4411K85.....\$3.11	4411K36.....\$3.70	4411K76.....\$2.20	4443SK16.....\$4.46
4411K26.....\$4.71	4411K86.....\$3.32	4411K37.....\$3.79	4411K77.....\$2.26	4443SK17.....\$4.53
4411K101.....\$7.04	4411K801.....\$4.32	4411K41.....\$5.88	4411K81.....\$3.20	4443SK21.....\$11.70
4411K102.....\$7.12	4411K802.....\$4.44	4411K42.....\$6.00	4411K82.....\$3.22	4443SK22.....\$12.55
4411K103.....\$8.16	4411K803.....\$4.80	4411K43.....\$6.08	4411K83.....\$3.40	4443SK23.....\$13.19
4411K104.....\$8.16	4411K804.....\$5.08	4411K44.....\$6.16	4411K84.....\$3.48	4443SK24.....\$14.09
4411K105.....\$8.28	4411K805.....\$5.28	4411K45.....\$6.16	4411K85.....\$3.48	4443SK25.....\$14.89
4411K106.....\$8.48	4411K806.....\$5.80	4411K46.....\$6.16	4411K86.....\$3.88	4443SK29.....\$14.89
4411K111.....\$9.60	4411K811.....\$5.89	4411K51.....\$8.18	4411K91.....\$4.36	
4411K112.....\$9.71	4411K812.....\$6.05	4411K52.....\$8.18	4411K92.....\$4.53	
4411K113.....\$10.25	4411K813.....\$6.55	4411K53.....\$8.29	4411K93.....\$4.64	
4411K114.....\$11.13	4411K814.....\$6.93	4411K54.....\$8.40	4411K94.....\$4.75	
4411K115.....\$11.29	4411K815.....\$7.20	4411K55.....\$8.40	4411K95.....\$4.75	
4411K116.....\$11.56	4411K816.....\$7.64	4411K56.....\$8.13	4411K96.....\$5.02	

SUBMICRONIC WOUND COTTON—0.5 MICRON RATING

Each	Cartridge Length	Each	Cartridge Length	Each	Cartridge Length
4411K981.....\$8.58	20"	4411K982.....\$11.80	30"	4411K983.....\$17.70	40"

Resin-Bonded Filter Cartridges

Resin-bonded filter cartridges are designed for high efficiency filtering of organic and petroleum oils. They are available in 10", 20", and 30" diameters. The ends are Double Open Ended (D.O.E.) unless otherwise indicated.

Micron Rating	10" Cartridge Length	20" Cartridge Length	30" Cartridge Length
150.....	4383SK18.....\$4.05	4383SK28.....\$7.58	4383SK38.....\$11.18
100.....	4383SK17.....\$4.33	4383SK27.....\$7.91	4383SK37.....\$11.40
75.....	4383SK16.....\$4.35	4383SK26.....\$7.96	4383SK36.....\$11.73
50.....	4383SK15.....\$4.35	4383SK25.....\$7.96	4383SK35.....\$11.73
25.....	4383SK14.....\$4.15	4383SK24.....\$7.04	4383SK34.....\$10.36
10.....	4383SK13.....\$4.24	4383SK23.....\$7.67	4383SK33.....\$11.78
5.....	4383SK12.....\$4.58	4383SK22.....\$8.60	4383SK32.....\$12.59
2.....	4383SK11.....\$4.97	4383SK21.....\$9.40	4383SK31.....\$14.22

McMASTER-CARR

Carbon Filter Cartridges

Spiral Wound Carbon

- Remove bad taste, odor, and sediment
- For point-of-use, single water faucet only
- Maximum temperature is 125° F • 4 1/8" cartridges have a 2 1/2" OD; 9 1/2" and 9 7/8" units have a 2 1/2" OD

Filled Carbon

- Spiral Wound—Choose from a standard cellulose media or a powder activated carbon media with metal adsorption resins for lead removal.

Cartridge Style	Micron Rating	Each
1 1/2" Cellulose.....	5	AS17SK15.....\$3.99
9 1/2" Cellulose.....	5	AS17SK4.....\$4.53
9 7/8" Activated Carbon.....	0.5	AS17SK6.....\$2.38
• Aneltek model CMR-10.....		

Pleated Filter Cartridges

Pleated Polyester

- Clean and reuse in many liquid applications
- Four square feet of media
- Maximum temperature is 140° F • 2 1/4" OD

High Volume

- Removes more solids than pleated units • Six square feet of media • Maximum temperature is 140° F • 2 1/4" OD

Micron Rating	9 1/2" Cartridge Length	20" Cartridge Length
50.....	4423SK13.....\$3.97	4423SK23.....\$7.38
20.....	4423SK15.....3.56	4423SK25.....7.38
5.....	4423SK17.....3.76	4423SK27.....7.38

Spun and Melt-Blown Polypropylene Filter Cartridges

Spun Polypropylene

- Removes sediment, sand, silt, scale and rust particles and has excellent chemical resistance
- Maximum temperature is 145° F • 2 1/2" OD

Melt-Blown Polypropylene

- FDA listed materials
- Use with food and beverages
- Absolute-rated 99.98% efficiency
- For consistent, reliable performance
- Maximum temperature is 200° F @ 15 psi • 2 1/2" OD

Micron Rating	10" Cartridge Length	20" Cartridge Length	30" Cartridge Length
40.....	4417SK19.....\$9.58	4417SK29.....\$17.57	4417SK39.....\$26.83
30.....	4417SK18.....9.58	4417SK28.....17.57	4417SK38.....26.83
20.....	4417SK17.....11.23	4417SK27.....20.70	4417SK37.....28.63
15.....	4417SK16.....14.53	4417SK26.....26.83	4417SK36.....32.57
10.....	4417SK15.....14.53	4417SK25.....26.83	4417SK35.....32.57
5.....	4417SK14.....21.76	4417SK24.....32.57	4417SK34.....32.57
1.....	4417SK13.....21.76	4417SK23.....32.57	4417SK33.....32.57

Stainless Steel Filter Cartridges

Standard

- Cleanable, reusable filters withstand most caustic fluids. These Type 304 stainless steel filters are fashioned by welding and crimping wire cloth; you'll never have to worry about the disintegration of sinter brazing or epoxy bonding. Filters are double open-ended with grommet seals that can't be accidentally dislodged. The maximum temperature for cartridges with Buna-N seals is 250° F. Cartridges with Teflon seals have a maximum temperature of 500° F. Maximum differential pressure is 80 psi; outer diameter is 2 1/4". Cartridges are 10" long.

Micron Rating	Filter Style	Buna-N Seals	Teflon Seals
100.....	Standard.....	4427SK81.....\$49.11	4427SK81.....\$76.29
100.....	Pleated.....	4427SK82.....51.54	4427SK82.....75.38
20.....	Pleated.....	4427SK83.....116.15	4427SK83.....140.00
5.....	Pleated.....	4427SK84.....187.89	4427SK84.....191.54

Inline Disposable Water Filter Cartridge

Remove unpleasant aftertastes in coffee, tea, citrus, and drinking water. Install on your water line—the activated carbon filter cleans up to 7/2-gallon of water per minute, leaving it practically free of sediment, unpleasant tastes, and odors. When the filter is used up, simply detach, discard, and add a new unit.

Available with 1/2" brass compression fittings for easy attachment to your 1/2" water supply line, or with a quick-disconnect fitting. Initial installation requires a quick-disconnect adapter sold separately. Color: Maximum pressure is 125 psi; temperature range is 32° to 125° F. Measures 2" Dia. x 12" Lp.

Filter with Compression Fitting.....	4420SK93.....\$10.18
Filter with Quick Disconnect Fitting.....	4420SK94.....9.17
Installation Adapter For 4420SK94.....	4420SK96.....6.65

1581