

MODEL 2502A TECHNICAL MANUAL

I. SPECIFICATIONS

II. OVERVIEW

A. System Considerations

B. Electrical Connections

III. INITIAL TESTS AND ADJUSTMENTS

V. TECHNICAL DESCRIPTION

A. Pulse Position Analyzer

6500 - 017

B. Interface Board

-[2^

C. Buffer Memory

-m .

II. SPECIFICATIONS

MODEL 2502 POSITION COMPUTER

Image Quality:

Resolution:	256 x 256 pixels
Linearity:	5% Relative displacement of any point in image. Output monotonic over entire image

Dead Time

300 ns per event
<400 ns dead time plus peak time, (effective dead time for paralyzable dead-time correction)

Background Rate:

<10 Events/s over image

Outputs:

Sum Monitor:
(On position analyzer unit)

Sum of all preamps for observing and calibrating MCP gain level. 0.4 Volts bipolar

X and Y position:

0 to 3.75 Volts analog
resolution of pulse position

Z axis:

Oscilloscope display blanking
(0 V on, 5 V blanking)

Strobe:

TTL pulse (300 ns) during each accepted event

Rate:

TTL signal that toggles with each input event exceeding "fast" threshold level

Digital Interface:

Digital interface directly compatible with 2420B PC/AT interface.

Power Requirements:

Selectable voltage ranges;
100/120/220/240 VAC, 50 Watts

II. OVERVIEW

A- System Considerations

The 2502A electronics unit receives charge pulses from a microchannel plate (MCP) and resistive anode position encoder assembly. The position data are stored digitally in an internal buffer memory. The memory contents are sent to the host PC/AT computer via a Model 2420B PCIO interface. Analog outputs are provided for monitoring pulse positions, MCP pulse amplitude and count-rate. A digital rate output is also provided. The preamp and position decoder electronics are packaged in a small box that should be mounted near the sensor. The power supply, monitoring circuits, FIFO interface and buffer memory are contained in a separate chassis unit. The typical system configuration is shown in Figure 1.

Vacuum feedthroughs. Most conventional feedthroughs can be used for the microchannel plate HV bias electrodes and resistive anode signal outputs. Although it is not necessary to use coaxial RF feedthroughs for the signal leads, it is important to minimize inter-electrode and stray capacitance between the signal leads.

Vacuum Lends. In systems operating at pressures greater than 10^{-10} torr that are not baked at temperatures higher than 150°C , low capacitance teflon insulated coaxial wire such as RG-187 A/U is recommended for signal and HV leads. Signal leads should be kept as short as possible and the total length from detector to preamp input should not exceed 15 inches. The signal lead lengths should not differ by more than 1 inch. Higher input capacitance will degrade resolution. In UHV systems, HV leads should be bare or ceramic insulated wire. Leads should be shielded from each other and from any other potential noise sources.

Grounding and Bypassing. The pulse amplifiers are extremely sensitive to high frequency signals, so careful grounding and shielding practices are essential. Other electronic equipment that is interconnected with the position computer electronics should be powered from the same electrical circuit, and should be isolated as much as possible from electric motors, SCRs and other sources of electrical noise. The signal lead shields should be grounded to the vacuum chamber at one point as near to the detector as possible. All components near the anode should be connected to this ground. The position computer electronics box should be electrically grounded to the vacuum chamber. One of the four mounting holes on the box bracket has the painted surface removed to provide a good electrical contact. The analog monitor output provides a convenient check on the input noise level. Some experimentation with different grounding techniques should be tried to find the lowest noise configuration.

Decoupling. The pre-amp inputs are direct coupled and cannot tolerate any DC voltage. If the resistive anode is operated at other than ground potential, it is necessary to decouple the anode with 1000 pf capacitors. (This isolation is standard on most SSL open face sensors.)

B. Electrical Connections

Sensor input connections. The leads from the vacuum chamber feedthroughs to the preamp inputs should be RG 187 A/U coaxial cable. These leads should be kept as short as practical and of equal length.

Analog monitor output. An analog monitor output is available at a BNC connector on the position analyzer unit. This output should be used only for test purposes and should not have a cable attached during normal operation. The sum monitor signal is the combined pulse amplitude from the four preamp outputs. This signal is used when adjusting the MCP high voltage and for measuring noise level.

Position analyzer to main chassis connections Two cables connect the position analyzer to the main chassis. Digital data are carried over the 37 pin connector cable, and the 9 pin connector cable provides the power interface.

Computer interface The computer interface is a 50 conductor cable that connects the main chassis to a 2420B interface board installed in the PC/AT computer.

Chassis signals Several output signals are available on BNC connectors on the chassis back panel. X and Y axis analog signals and a Z axis strobe signal can be connected to an oscilloscope to provide a real time display of the image. The logarithmic count rate and MCP level signals that can be selected for video display are also available on back panel connectors. Digital event pulses (TTL level) representing both total and accepted events are provided.

III. INITIAL SETUP AND ADJUSTMENTS

With the detector installed in a vacuum chamber, use the following procedure for the initial turn-on. The chamber should have a signal source such as an ultraviolet lamp, ionization gauge, or filament (a flashlight bulb filament makes a convenient source). It need not be evacuated for the preliminary checks.

1. The detector high voltage should not be turned on in pressures higher than 10^{-5} torr.
2. With an oscilloscope connected to the "sum monitor" of the position electronics unit, turn on its power (detector HV remains off). A few millivolts of "white" noise should be visible. Check for any spikes and transients. The threshold for accepting events is 60 mV, so all noise should be substantially below that level. Try to improve grounding configuration if noise is excessive.
3. Evacuate chamber and slowly turn on the detector HV supply to a convenient low voltage (100 - 500 V). Check for any additional noise on the "sum" monitor. If the HV supply is properly filtered and isolated, no change should be seen. If any large increase in noise level is observed, check for possible sparking or discharge in the system.
4. Slowly increase the HV supply to 1700 V while continually monitoring the "sum" signal. Stop if there is any sudden jump in noise level. Once the HV is set, illuminate the detector with a small flux of charge particles or photons. Pulses should be visible above the noise level on the "sum" monitor. Signal pulses have a bipolar shape with a total pulse width of about 250 ns. Adjust the HV until the mean pulse height is about 200 mV. The pulse electronics will clip pulses at about 400 mV. It should not be necessary to increase the HV above 2 kV for this adjustment. The "input level" meter should read near mid-scale when the MCP gain is correct.
5. Check the ^Hposition monitor¹ to see that pulses are being detected. The position output should be jumping to levels between 0 and 3.75
6. Turn off the illumination source and note the number of background counts (i.e. "sum" pulses >0.06 V). This rate should be 10/second or fewer. If the rate is substantially higher than this, there may be a source of electrons, ions or photons in the chamber. Possible spurious sources are corona from a HV lead, ion gauges, or a source of ionizing radiation. Observe the spatial distribution of the background counts on the "position monitor"^H (a storage oscilloscope is

convenient for storing these low rates). These events should be randomly distributed. A concentration of counts in one area indicates either a "hot spot" on the MCP, or electronic transients being picked up by one or both preamps (a common manifestation of transient pick-up problem is an apparent "hot spot" in the center of the image due to a signal picked up by both preamps). A test that usually distinguishes transients from real hot spots is to decrease the high voltage of real MCP pulses to just below the discriminator threshold and disappear.

IV. TECHNIQUE OF THE HSCPHION

A. Pulse Position Analyzer

The Model 2502 detector system includes a resistive anode image converter, shaper amplifiers and position analysis electronics, and a buffer memory with computer interface. The resistive anode is placed behind two microchannel plates having approximately 5×10^6 gain, so that each incident particle produces a charge pulse on the anode that can be analyzed by the position electronics. The operational function of each component is described briefly in the next section and is followed by a technical description.

1. The resistive anode is a diffusive RC transmission line that divides charge between output electrodes in proportion to the proximity of the input pulse.
2. Charge sensitive preamps and shaper amplifiers convert the low level charge pulses into high level bipolar signals suitable for driving the position analyzer circuit. The amplifier shaping network has been matched to the resistive anode for optimum signal to noise ratio, short pulse duration, and low baseline shift.
3. The position analyzer electronics generate a digital output proportional to the pulse position coordinates of the form

$$X = B/(A + B)$$

The position analyzer logic selects acceptable pulses that have amplitudes between the upper and lower level discriminators. Accepted events are stored in the output register and an output strobe signal is generated.

Circuit Details

The major functional blocks of the pulse position analyzer are the four input amplifiers, a threshold amplifier, summing power amplifiers, ratioing A/D converters, and control logic. The following description refers to the schematic diagram.

The input amplifiers consist of transistors Q1-Q14, Q19-Q32 and their associated resistors and capacitors. The four amplifiers are identical so this description refers only to the "A" amplifier. The input stage is a common base charge amplifier followed by a two stage emitter follower. The second stage is a four transistor voltage amplifier (Q5-Q8) that drives the main LCR shaping network. A second voltage amplifier (Q9, Q10) drives the second shaping network. The final stage is a voltage amplifier/driver (Q11-Q14).

A Hep... amplifier with very C... swiping $i_{uhc} < 1$ for Lho thrcKholej signal. The output of the four channels is amplified and amplified by Q15-Q18. The output pulse of this amplifier has a duration of about 50 ns, which is much shorter than the main amplifier pulse duration and facilitates effective pulse pile-up rejection by the control logic. The threshold level is set by the comparator U6, which also comprises a one-shot multivibrator with a time constant of 50 ns. Two identical power amplifiers are used to drive the inputs to each of the ratio A/D converters U9. The "numerator" or "A" power amplifier is made up of Q33, Q34, U1, and half of U2. The transistors and the LH0002 comprise a unity gain pulse amplifier having low output impedance. The LM13600 provides baseline restoration to maintain a stable de output level for the pulse amplifier. Diodes D15 and D16 clamp the positive and negative output swing to protect the ADC. The "denominator" signal is formed by summing "ABCD". The power stage for this signal consists of Q35, Q36, U3, and half of U2.

The actual position conversion is carried out by the A/D converter, U9. This is a flash ADC containing 256 separate voltage comparators attached to taps on a voltage divider resistor string. The analog input signal is compared with this reference voltage divider, and all comparators at taps below the input voltage level are triggered. The comparator outputs are coded into a standard binary format. The digital output represents the ratio of the analog input and reference input signals. The ADC must be strobed when the analog signal pulses are near their maximum value and the timing for this strobe signal is generated by the control logic.

The control logic selects events for analysis and generates the required timing and strobe signals to operate the ADC and output register. To be accepted, a pulse must have an amplitude lying between the threshold and overload levels and must not overlap other pulses. Logic pulses from the "threshold" comparator trigger the "peak delay" and "dead time" one-shots. The end of "peak delay" clocks UIOb, which starts a "convert" pulse and also triggers the "convert" one-shot U11a. This first convert pulse actually initiates the A to D conversion by U9, but a second convert pulse is needed to transfer the data to the output buss of U9. The "convert" one-shot generates this double pulse by a combination of its trigger delay and the normal timeout. After a delay of about 50 ns the "convert" one-shot is set, inhibiting the "convert" signal. After an additional 40 ns the one-shot times out and the inhibit is removed, producing the second "convert" pulse.


If a second threshold pulse arrives before the end of "peak delay" UIOb will be set and no "convert" pulses will be produced when UIOb is clocked. Once the A/D conversion has been initiated the output strobe will be generated if no "overload" condition is detected. The overload comparator is polled by U14 about 20 ns after the end of "peak delay". If "overload" is false, U14 will be set and subsequently reset by the timeout of "dead time". The falling edge of U14 triggers the "output strobe" one-shot U11b.

B. Interface board


— /x progress (6500 - a1?) 2

C. Buffer memory (see 2412 buffer memory manual)

(6500-0189)



2502 POSITION COMPUTER
 BOARD A VERSION S2D
 SURFACE SCIENCE INSTRUMENTS
 C.M. CARLSON 10/7/83



2502 POSITION COMPUTER
 BOARD B S/N 9337
 SURFACE SCIENCE INSTRUMENTS
 C. M. CARLSON 10/77/83

DIODES:
 1N914A

UNMARKED TRANSISTORS:
 PNP 2N3904 NPN 2N3904

