Instructions for Hawk Software Service Physics, Inc



March 2010

Highlights of Changes form ESCA 200X to ESCA Hawk V1.01

Spectrum Display Window (Graphic Viewer)

The Spectrum Display Window is used in both the capture and analysis programs. The following additions are included:

- Create stored preferences for appearance of spectra.
- Individual appearance settings for
 - o Graphic Viewer
 - o Clipboard
 - Printing only available in Analysis application
- Quick Titles and Figure descriptions can be added to Spectra.
- Icon to save Spectra to file
- Subtract back ground/Restore back ground icons
- Icon for adjusting x,y scales

Capture Application

These features are in addition to the Spectrum Display Window.

The following additions are included:

- Direct transfer of parameters between the MRS table and the ESCA control Panel
 - Test Region button transfers the parameters of the marked region of the MRS table to the ECP and starts the data capture.
 - Update Region transfers the parameters in the ECP to the marked Region in the MRS table
 - A blank Region in the MRS table can be filled. Automatic assignment to appropriate MRS table functions is performed.
- Elapse time estimate provided for MRS table.
 - Scan time is provided for each region
 - o Total time for all regions is provided
 - Time to run from marked region to end of MRS table is provided.
 - During Capture the total time for the current region to the end of the MRS Table is provided.
- MRS Table improvements
 - Each region Time/Step may be programmed in 25 millisecond steps.
 - In UnScanned mode, Detector Width and Channel Width are displayed.
- ESCA Control Panel (ECP) Improvements
 - When ECP is closed the state of the parameter set is stored as the "Last Used".
 - At any time you can store the current state of the parameter set as "Default".
 - You can set the ECP to open using the "Default" or "Last Used" parameters.

- A Count Rate mode is added using an exponential decay with time. This provides monitoring the effects of changes in charge compensation, crystal adjustments and spectrometer adjustments with and adjustable time filter.
- "Quick Shift" Binding Energy change. This control behaves like the "brush" mode on IPod displays. The control will move the CBE as much as half a window with in either direction.
- Toolbar icons for quick opening of "Analysis", "ESCA Control Panel", "Big Display" and "Export" of spectra.
- Setup ESCA (Found in the Settings menu)
 - Configuration table Added buttons to allow saving and loading the configuration table. Use to back up calibration information.
 - Added a default peak fit to locate peaks in Detector Width Calibration
- Add additional diagnostics to detect connection problems with detector.
- Improved reliability of Depth Profile setup, running and ion gun control.
- Region Name cell in MRS table is used in Quick Titles for auto assigning names.

Analysis Program

These features work with the new features of the Spectrum Display Window

- Fast Report (FR), a powerful report generation program, has been imbedded for hard copy output of analysis information.
 - A set of templates are included to provide standard print configurations.
 - The full Fast Report Designer is imbedded to allow modification and storage of templates.
 - Fast Report (FR) provides extensive tools for adding spectra and data to predefined reports.
 - Report Review provides options for a wide range of export of assembled data.
 - Tables can be automatically added to Fast Report
 - Column Titles, Placement, Visibility, Format and Font are selectable in Fast Report.
- Editing of Project, Experiment and Region names.
- Delete for Projects, Experiments and Regions.
- Print Tool Bar
 - o Select Template
 - Add (one) Region to Fast Report (+Regn to FR)
 - Add all Regions to Report (+ All Regn to FR)
 - Empty Fast Report (Empty FR)
 - Review the assembled report (Review Report)
 - Print out assembled report (Print Now)
 - Design or modify a template (Design Template)

- Table Functions
 - Export Tables to EXCEL
 - o Copy full table Clipboard
 - Copy single row or column to Clipboard
 - o Launch Table with cells unprotected

DP Viewer

DP Viewer can be launched directly from the Analysis Application tool bar.

- Tool Bar
 - Add Control of line color
 - Export all cycles of one Region as EMF files. Each EMF file contains one spectrum.
 - o Zoom Control
 - Support for Predefined Chart properties after export to EXCEL

SEE - OVERVIEW OF NEW FEATURES after next section for examples of these features.

Recommended system configuration.

- 1. Pentium IV or higher, 256K Ram, 1Ghz
- 2. USB 2.0
- 3. Windows XP operating system with service pack 2
- 4. Microsoft Office 97 or higher (Word, Excel and Power Point).
- 5. 19" Flat panel Display.
- 6. CD R/W
- 7. Color Printer

Windows 2000 operating system **MUST HAVE SERVICE PACK 4 INSTALLED** to upgrade to ESCA Hawk. The export to Excel requires various security configuration changes in Excel to allow the macro program to operate. If problems occur exporting to Excel see notes at end of install section or contact Service Physics for help.

Preparation of system for upgrade.

- Backup calibration file. Older ESCA software does not have a calibration export utility. The ESCA Hawk CD has included a program that will create a backup calibration file. The utility is in the Folder, "Hawk Support" > "ESCA Calibration Config". Copy the file, "XpsRegReadSetup.msi" to the existing ESCA Program folder where the old ESCA program resides. Double click the file to start the utility. Follow the directions to create a backup of the calibration configuration. We suggest you use a name such as ESCA Cal 080521. This provides the date code year 08, month 05, day 21. ESCA 25 and ESCA Hawk have this utility built in. You will be able to load the calibration file after ESCA Hawk is installed. A history of your calibrations can be created from the backup files saved in the Configuration subfolder of you ESCA Hawk program folder. See page 22 in this manual.
- 2. Uninstall existing ESCA program. Do not uninstall MDAC, National Instruments Software, ACCESSRT, CH Products or Motion Planner. See section below if you are updating your operating system.
- 3. Updating Operating systems. If updating system from NT the IEEE488 National Instruments card will not work. A new PCI GPIB card needs to be purchased. The IEEE488 version 2.22 (or higher) driver should be used. If updating from Win 2000 to XP it is best to reformat the Hard Drive and install a new copy of windows XP with service pack 2. The PCI GPIB card used on 2000 will work but the Version 2.22 (or higher) driver must be used and not the old driver that may be on the National

Instruments CD. IEEE 488.2 V2.24 is included on the ESCA Hawk CD. Updating requires additional memory and will run much better with at least 1 GHZ processors.

- 4. ESCA 2000 V1.02.0x releases do not have the backup utility discussed in paragraph 1. It is very important to follow the backup instruction if your are changing operating systems or using a new computer. If you are making no changes to the computer operating system the old configuration should not be lost. It is still advised that you follow step 1 in case of computer problems.
- 5. There have been problems uninstalling ESCA2000 V1.02.0x. A utility is included on the CD called Uninstall Cleanup. This is found in the Uninstall Utilities folder. This program cleans up a number of problems with the registry and removes some components that are not removed by uninstall. This program should be used to clean up any installation that has had previous problems including ESCA Hawk V1.01. Uninstall Cleanup provides a report about the cleanup. The desired result is that each of the components is not found. If found it uninstalls them and deletes the file from the system 32 folder. Please have the utility make a log file if any cleanup is needed. Please e-mail the log file to us.

Install upgrade.

- 1. If installing ESCA Hawk on a new computer or on a computer with newly upgraded operating system then do the pre-install in step 2 below. If no change to the operating system then jump to step 3 below.
- 2. Pre-install support components on new computers and computers with updated operating systems. See table 1 below for components that need to be installed.
- 3. Install the new ESCA Hawk program. Put the new CD in the disk drive. Find the ESCA Hawk folder and run the Setup program. Accept all default options. Accept keeping newer components that may be on the hard drive.
- Copy the Calibration Configuration file, saved in step 1 of the "Preparation of system for upgrade" section above, to the Calibration folder in the ESCA Hawk directory.
- 5. If this is a new installation you will need to get a Serial Number form Service Physics by e-mailing either <u>bob@sphysics.com</u> or <u>barb@sphysics.com</u>.
- 6. Open the Capture program. Select the settings menu. Select Setup ESCA. Open the Configuration Tab and select the Load from File button. Navigate to the Calibration folder and load your backup calibration parameters.

		OT LET I DI LD	
VERSION	DATA BASE	MOTION	JOY STICK
ESCA HAWK A	ACCESSRT*	NONE	NONE
ESCA HAWK B, BU	ACCESSRT*	NONE	NONE
ESCA HAWK G, GU	ACCESSRT*	NI 488.2 V2.24	NONE
ESCA HAWK E, EU	ACCESSRT*	6K4 MOTION	CH PRODUCTS
		PLANNER	

TABLE 1 - SUPPORT COMPONENTS FOR ESCA HAWK.

NOTE: USB VERSIONS WILL AUTOMATICLY INSTALL THE USB DRIVER. LEAVE THE USB CABLE DISCONNECTED DURING THE INSTALLATION OF THE ESCA HAWK PROGRAM. AFTER THE PROGRAM IS INSTALLED THEN CONNECT THE USB CABLE. THE WIZARD WILL ASK IF YOU WISH TO LET WINDOWS FIND THE DRIVER. BE SURE IT IS NOT ASKING TO USE THE INTERNET. IT SHOULD LOOK ON THE CD FOR THE DRIVER. WINDOWS WILL AUTOMATICLY FIND THE DRIVER ON THE CD AND INSTALL. IF THE USB DRIVER IS ALREADY INSTALLED WINDOWS WILL SKIP THIS STEP.

*ACCESSRT IS FOUND IN THE SUPPORT FOLDER ON THE ESCA HAWK CD.

OVERVIEW OF NEW FEATURES

New MRS Table design

		Update Region	Tesi Regio	t on	• (•)	Over Write Append	Run ⊂ One + End _⊂ =	0:02:19	3	Full Run Row to E	Time Ind Tm				
	1	Function			LBE	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time		-	Clear
		Survey		-	0	250	800	1	4	1	100	0:00:34	1		
		G Survey	1000	-	1										1
	2	Function			CBE	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time			
		ResHi		•	284	20	300	1	2	0.065	100	0:00:48			•
		G HR R2		•	1										•*
	3	Function			LBE	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time			A
		HiSen		•	50	100	800	1	4	.4	100	0:00:38			X
		G Sen R4	1	•	1										X
	4	Function			CBE	Det Width	Spot Size	Time	Res#	Ev/Chnl	Tm/Step	Scan Time			Ŧ
		UnScann	ied	•	84	21.500	800	00:15	4	0.168	100	0:00:19			
		GUS		Ŧ	1										Ŧ
		END T													
1													•		

- Display Detector Width for UnScanned Function
- Display Channel Width of the Detector for UnScanned Function
- Time/Step (Dwell Time variable) may be input in 25 millisecond steps
- Actual elapsed time for scan. This is longer than the time spent capturing data due to computation overhead. Longer Dwell times reduce percent overhead.
- Elapse time for complete MRS "Full Run Time"
- Elapse time from current active row to end of MRS "Row to End Tm"
- "Row to End Tm" tracks the progress of a running MRS table.

Test Region



• Selecting "Test Region", in controls above the MRS table, transfers data from active row of the MRS table (Yellow Highlight button) to the ESCA Control Panel. If panel is not open it will open, down load parameters and start data collection.

Update Region



- Changes to ESCA Control Panel are uploaded to MRS table.
- Resolution was changed from 2 to 1, Number of scans from 1 to 2 and the Scan Time was updated.

Update Region used to build MRS table

Spectrometer Controles	- Initial Param	1	1	Function	L	.BE I	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time	-	Charles 1
Aperture	lastlised	Ш		Survey 🗖	- 0) ;	250	800	1	4	1	100	0:00:34		
Flood Ion X-ray Gun None V	Last Coold	Ш		G Survey 1000 🗖	- 11	0									•
dun dun operate i	C Default	Ш	2	Function	C	BE \	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time		
	Cat Commit	Ш		ResHi 🗖	- 2	284 3	20	300	2	1	0.065	100	0:01:31		•
CBE Spot Size Res# Cap I	as Default	Ш		G HR R2 🗖	- 1	1						_			•*
204 200 1 1		Ш	3	Function	L	BE \	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time		A
	Scans	Ш		HiSen	- 5	50	100	800	1	4	.4	100	0:00:38		×
Scans Det Width Ev/Chol Tm/Sten		Ш		G Sen R4 🗖	- 11	0									X
		Ш	4	Function	C	CBE	Det Width	Spot Size	Time	Res#	Ev/Chnl	Tm/Step	Scan Time		不
2 3.641179 2.844 100 💌		Ш		UnScanned 📃	- 8	34 ;	21.500	800	00:15	4	0.168	100	0:00:19		
		Ш		G US 🗖	- 11	1									_ ₹_
	Count Rate	Ш	5	Function											
STABT * Abort Exit		Ш		•	-										
		Ш		<u> </u>		1									
										END				-	
			•												

Add a blank line to MRS table.

Spectrometer Controles	Initial Param	1	1	Function		LBE	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time	-	Clear
Fland In You Cur Aperture	C Last Used			Survey	Ŧ	0	250	800	1	4	1	100	0:00:34		
Gun Gun Operate None -				G Survey 1000	Ŧ	1									1
	C Default		2	Function		CBE	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time		
				ResHi	Ŧ	284	20	300	2	1	0.065	100	0:01:31		•
CBE Spot Size Res# Cap <u>T</u>	Set Current			G HR R2	•	1									•
			3	Function		LBE	Wnd Width	Spot Size	Scans	Res#	Ev/Step	Tm/Step	Scan Time		A
	Scane			HiSen	Ŧ	50	100	800	1	4	.4	100	0:00:38		X
Come Det) (idth Eu/Chol To/Stop	Juans			G Sen R4	•	1									X
1 Scans Det width Ewichni Thivstep			4	Function		CBE	Det Width	Spot Size	Time	Res#	Ev/Chnl	Tm/Step	Scan Time		T
2 3.641179 2.844 100 -				UnScanned	•	284	3.641	300	15	1	2.844	100	0:00:19		
, , , , _				GUS	•	I									. ★
,	Count Rate		5	Function		CBE	Det Width	Spot Size	Time	Res#	Ev/Chnl	Tm/Step	Scan Time		
				UnScanned	Ŧ	284	3.641	300	15	1	2.844	100	0:00:19		
				GUS	•	I									
										END					
														•	
			4											<u>ا ا</u>	

Selecting Update Region creates the correct function and enters the parameters from the ESCA Control Panel.

Time average with exponential decay.





ESCA Control Panel Initialization

ESCA Control Panel	
Spectrometer Controles Flood I Ion X-ray Gun None I Gun Gun Operate	Initial Param Initial Param Last Used Default
CBE Spot Size Res# Avg I 89 300 ✓ 4 ✓ Scans Det Width Ev/Chnl Tm/Step	Set Current as Default Scans
2 21.5 0.167 100 -	562
START * Abort Exit	5521595

- Set Current As Default Stores the current parameter set.
- Each time the ECP is closed the current parameter set is saved as "Last Used"
- The selection buttons allow the user to chose which parameter set will be used the next time the ECP is opened.
- The parameter set is the eight parameters in the black box starting with the CBE and ending with the Tm/Step.

CBE/LBE and Avg T/Cap T Buttons

The buttons above the start BE box is fixed to CBE for UnScanned spectra. You can select CBE or LBE for scanned spectra.

The button over the Time box can be switched from Cap T (Capture Time) to Avg T (Averaging Time) in the Unscanned Mode. It displays the elapse time for a scanned spectrum.



The Project, Experiment and Region information can be edited using the MDB Editor. The Editor can also be used to delete at any level of the project tree. The Data Base structure is very detailed and includes reduced data sets that have been stored. Please keep in mind that deleting an Experiment or Project deletes all information stored for the project or experiment.

Edit Names and Descriptions		
Project:a Project:A Quickie Project:D Project:Test B-Raw MRS Regn4(73.3.94.8 eV, 128 pts) Regn3(50.150.4 eV, 251 pts) Regn2(274.294.1 eV, 309 pts) Regn1(-251. eV, 251 pts)	Edit Experiment Information Project:Test Experiment:4 functions Has 4 spectra and 0 data reductions. Change project of this Experiment (moves to or creates a new project): Test Change Experiment Name To: 4 functions Change Experiment Description To: short scan Change Operator Name To: Bob Delete Experiment	
		Quit

Edit Names and Descriptions		
 Project:a Project:A Quickie Project:D Project:Test Experiment: 4 functions Braw MRS Regn4(73.3:94.8 eV, 128 pts) Regn3(50150.4 eV, 251 pts) Regn2(274294.1 eV, 309 pts) Scan on dark spot(251. eV, 251 pts) 	Edit Region Name Region: Scan on dark spot SpectrumName: Scan on dark spot LBE: 0; WW: 251; Npts: 251 Scans: 1; Cycle: 1 Operator: Bob Comment: short scan Change Region Name: Scan on dark spot	
	Delete Region Update Region	
		Quit

Edit Names and Descriptions		
 Project a Project A Quickie Project D Project T est Experiment: 4 functions Regn4(73.3-94.8 eV, 128 pts) Regn3(50.150.4 eV, 251 pts) Regn2(274294.1 eV, 309 pts) Scan on dark spot(251. eV, 251 pts) 	Edit Project Name Summary for Project: Test Contains 1 experiments, with 0 data reductions, and with 4 spectra. Change Project Name To: Test Delete or Update removes or updates the entire Project and all Experiments and Spectra. They can not be recovered. Expand Project tree and select Experiment to modify the Experiment Name or Experiment Description. Seperate Operators names can be input for each experiment. Individual Region names can be updated by selecting the Region in the Raw or Analyzed data sets. Delete All Update All	Å
		Quit

Spectrum Display Window (Graphic Viewer)



Create Spectra Appearance Preferences using the Chart Preferences Dialog

Chart Preferences Dialog



There are three targets window types; workspace, print and copy to clipboard. You can create a preference configuration using the "Create New…" button. You can edit any configuration using the "Edit Chart Preference Editor" button. One configuration must be checked for each target.

OPEN Chart Preferences Editor

Lachart Preferences Editor	
Chart Perspective	
Grid	2,400
Visible Advanced Editor	2,200
Line Style Panel Color	2,000
C Solid	1,800
© Dot	(V 1 400
C Dash Dot Dash Dot Dot	1,200
C Small Dots	Х 1,000
	800
	600
Horiz axis label BINDING ENERGY - EV	400-
Vert axis label COUNTS	200
	599.6 549.6 499.6 449.6 399.6 349.6 299.6 249.6 199.6 149.6 99.6 49.6
	BINDING ENERGY - EV
Cancel Done	

Simple changes of spectrum are selected and reviewed on Chart page.



Change grid type and visibility. Change the background color. Change text label for x and y axes.

Note: The line color is changed using the color wheel icon on Graphic Viewer Toolbar. Line colors are unique to the spectrum not the window preference.

NOTE: WE HAVE INCLUDED THE NEXT EXERCISE TO INTRODUCE THE ADVANCED EDITOR.

PLEASE DO THE FOLLOWING EXERCISE TO SET UP YOUR NEW SOFTWARE!

Exercise 1. ADJUSTING THE CLEAR SPACE ABOVE THE HIGHEST PEAK. Spectrum with highest peak touching top of window.



Advanced Editor



The Axis tab is selected.

To change the clear space at top of spectrum select:

Advanced Editor >Axis >Position. Change the value in the "Start %" entry box.

When the Advanced Editor is used the settings will be saved when "Done" is selected on the Chart Preferences Dialog. This setting will be used for the target window (Workspace, Print or Clipboard) selected. IF YOU WANT THE SETTING TO MATCH IN ALL TARGET WINDOWS YOU WILL NEED TO EDIT EACH CONFIGURATION.



Quick Titles and Figure descriptions

Quick Title Entry for Spectrum dialog



Result of using Apply All Button

Region name is established one of 4 ways.

- By default when spectrum is collected. It is the line number in the MRS table.
- By entering a name in the region cell of the MRS table. See Below.
- From the Region Name entered in the Capture, Depth Profile table.
- By editing the Region Name in the Scan Table.
- The Tools for Smooth, Differentiate, Integrate, Add and Subtract provide entry of Region names for the derived spectra.

Examples of Region Name entry



Region name entered in MRS table in Capture Application is used in Analysis Scan Table

Depth profile Region Name

	5	De	epth Profile	Dialog							
Г	Template Name					e Date	Last Modifie	-	61		
Au on SS				7/21/	/2000						
						_					
2	òΡ	ot S	ize Etc	n	Flood	Gun	Cycles				
	40	0	00	00:15	off		6				
l					Γ					1	
F	lo	tate							Set As	Default	
ſ	ofl	÷						_	Yes		
l											
L					Region	Definitio	ons				
L			RegionName	CBE	Scanned?	ww	Scans/Time	Res#	EV/Ster	- 💾	
L		►	Au 4f	84		20.92	5	4			
ł			C1s	284	V	20	1	4	0.4	- 🏹	
	Enter Region Name as peak identification. Links to Auto Scofield lookup in DP Viewer.										
				2 <u>C</u> ar	ncel <u>D</u> ele	te <u>F</u>	aste				

Edit Region Name in Scan List



Region can be renamed in scan table.



Tools auto create new Region Names.

Area	Peak Au	to Eind	🎗 Dp 🤆	🔪 🛛 Adj Axis 🗌 🛓	⊾ ≁ ∉	🍯 🖻 🕹	<u>k</u> [۵ 🗵 💽	Title	Proj Name: 1Agilent	
OnMark	150(Smooth)	:								Exp Name: NW Ribbon 0	FF Marks A
					-		-			Exp Desc:	
					1		1			Analysis: Raw data	
				1	1		1			Scan	\ \
	S .	Quick Title Ent	ry for Spectru	m						OnMark150(Smooth)
220 -										Off Mark 150	
200-		Title will be user	for all regions of	data set						UnMark150	
400							_				
180-								Add Title	Clear Title		
¹⁶⁰ م		1									
노 140 -		🔽 Auto Include	e Region Name	Exam	ple - Reg1:	Your Title					
ର ମ ସହ											
100		Figure Descripti	on for this Spectra	1					/		
							_				
80-								Add Foot			
60 -		Example	Cleant	a + 1 - a			4:41.	and f		Leone Auto	
40-			Clear t	ext boxe	es lo I	emove	uue	e and ro	boler lext.	Leave Auto	
			Inc.che	ecked I	Regio	n is au	oma	ntically	used as he	eader	
		Apply A	ine enc	encu. I	ce gio	11 15 uu	.01110	urearry	ubeu ub me		
Idle	eV-27	7.88 Coupts: 46	CL2+ 0+ 4d5								
Idle	eV:27	7.88 , Counts: 46	U12s Os 4d5								

Region as only title.



Modify Header and Footer fonts

Select: Chart Preferences > Chart Preference Editor > Titles > Fonts > Size and Color.



Region header with Font, size and color changed.

See section on Fast Reports (Page 23 - 32) for examples of using the Graphic Display Header and Footer information to enhance the Fast Reports.



Notice the older font, size and color are in use. The graph background was set to white and the grids are present. To bring them up to date we will use the Chart Preference Dialog to edit the Default Chart Copy preferences. After preferences are set they will remain even if the Application is closed. They are fixed until you change them. If preferences are entered in the Analysis Program they will be installed in the Capture Program the next time it is opened.

Icon to save spectrum to File



Customize Background Color for Power Point Presentations

Preferences for clipboard are used for Export. Clipboard preferences set for pasting into Power Point.



Clipboard was pasted into this instruction document in the above document. Now each Region selected will have the same format. Region names will be automatically entered as the title.

Demonstrating line color control and auto assignment of Region name as title.





Background removed/Background restored.

Background remove and restore icons.

Axis Adjusting Tool



If the Graphic Viewer Window is small there will not be enough room to place the lettering in the auto fine format. Auto Fine and Auto Coarse have different spacing models but in a small window they default back to the Graphic Viewer model. The Fixed Model lets the user set the increment between labels. These fixed values will be used regardless of the Maximum and Minimum of the data. The fixed scale model will produce "bad behavior" during data collection.

Capture Application

Settings Menu

🛢 Esca Capture (C:\Data Fol	ders\ESCA 25	oVAgilent.	mdb)	
Eile View Settings Help Exit Image: Settings Set System Regist Capture Settings Image: Settings Edit Scofield Table Image: Settings Edit Scofield Table Image: Settings Set up ESCA Image: Settings Create a new Sign Image: Settings Capture Settings Image: Settings Set up ESCA Image: Settings Create a new Sign Image: Settings Default Image: Settings DP Tess	ry Defaults hature file	o Set up	SPOL	Ttility
📓 XProbe Spectrometer Setup				
Spectrometer Control Panel	Acquisition Configuration	DetectorWidth	DAC V1 Curves	
Start Accumulation CenterBe 278	Registre mes controlling ages will not become) configuration are dis permanent until the c	splayed below. You ma ommit button has beer	ay edit them but n pressed.
WindowWin Configuration tab	Refresh	Commit	Save to File	Load from File
Nscans Nrea	Entry	Value	Type	Notes
Resolution	InterruptNumber	10	DWORD ISA int	errupt number used by buffer me
HorizLinearization Dithered Signature	GPI012Base GPI013Base GPI014Base	NEW	- Save and	l load
Play Mode	Spot Type Ctl12Null	Canora		Suration table.
BindingEnergy 278	Ctl13Null	To		30
Resolution 1	Ctl14Null	16	DWORD Control	s logic of strobe handshaking bi

"Save to File" and "Load from File"

The calibration information in the Configuration Tables is extracted to a XML file to provide backup of the instrument calibration. The "Load from File" button loads the backup data into the operating system registry to restore you system calibration. Older software did not have this capability.

The CD provided has a folder labeled System Calibration Utility. This program should be loaded on older installations and run to provide the backup file. After the upgrade the "Load from File" utility can be used if the calibration information becomes lost. This must be done if a new operating system is used in for the upgrade.

The "Save to File" button should be use to save any changes after calibration or service to the instrument. The file should be saved to provide a calibration history.

Hawk Analysis Application

Fast Report

Fast Report is a robust, imbedded program that connects the Graphics Viewer, Data Spread Sheet and Capture Parameter list to a predefined template. The assembled information on the template can then be printed, stored or exported as HTML, EXCEL, RTF or PTF files. The template designer provides a flexible way to customize your output. Predefined text can be incorporated in the template to provide standard reports with the spectra and graphs automatically included.

We have included 7 standard templates to provide convenient printing and exporting. To print out a spectrum you only need to select "Print Now" and the active spectrum will print. To print a spectrum with the composition table choose the appropriate template and select "Print Now". If you have a large number of spectra in the region table you select "+ All Regn to FR" (add all regions to Fast Report) and then select Print Now.

You can always preview the print out before it is printed or exported. It is not possible to edit the information once it is entered into Fast Report. It is important to have the spectrum and tables the way you want them before adding to fast reports.

One spectrum, with or without a table, can be added to Fast Report at a time by selecting "+ Regn to FR".

The report template can be Empted (Cleared) at any time.

The following examples demonstrate the flexibility of this printing system.

Selecting a Template



Dialog used to select template. Templates are available that show Spectra with header, collection parameter list, and data tables. The header provides a space for Laboratory name and logo. Custom templates can be created, stored and retrieved using the Browse button.

The following pages display examples or printing with the seven templates.





CAUTION: THE COMPOSITION TABLE COMBINES DATA FROM ALL SPECTRA. THIS IS TO ALLOW MIXING OF HIGH SENSITIVITY SCANS TAKEN AT RES 4 WITH SURVEYS TAKEM AT RES 4. IF TWO SURVEYS ARE IN THE SCAN TABLE BE SURE ONLY ONE TABLE HAS "Y" CHARACTERS IN THE INCLUDE COLUMN. YOU MUST ALSO BE SURE YOU MATCH THE SPECTRA ADDED TO FAST REPORT WITH DATA TABLE INCLUDED IN FAST REPORT. THERE IS ONLY ONE COMPOSITION TABLE AND THE "+REGN TO FR" BUTTON WILL ADD IT TO ANY ONE OF THE SURVEY SPECTRA IN THE TABLE.

PEAK FIT TABLES ARE MATCHED WITH THE SPECTRA USED TO COMPUTE THE FITS. The above caution does not apply to the peak fit tables. This can be seen in the example above.

The "Quick Title" tool button allows information to be added to the fixed Templates. See pages 12-16 for details on setting up Quick Titles.

Special print-to-file options are provided in "Review Report" utility. The graphic can be printed, sent to the clipboard, or exported in a number of formats.

Table Functions Tool Bar



The clip board icon picks up the table. In this example the table was pasted into this word document. The table was then edited in word to eliminate the sensitivity column. Words table editing tools provide column, row and cell formatting and editing.

XPS Line	Adj'ed Be	Norm Area	Atom %
O 1s	530.062	366.010	53.830
C 1s	284.556	154.538	22.728
Zr 3d	181.965	159.395	23.442



The tables can be exported to EXCEL or the embedded spreadsheet can be used for copy and paste operations. The "copy selection" icon can be used to pick up a single row or column. You can select multiple rows or columns if they are continuous without skipping any rows or columns.

	tale Recetor Administrator Date 9/20/2002									
Load Data L Add Data										
	INumber Prints in Q ITemplate Name Analyst Name									
je comp:POS1	Eile • 😰 Shift 12 📴 Copy Selection 👪 Launch Spreadsheet									
je comp:PUS2 Zr02	AA4									
ZrO2 Multi Region	1 XPS Line Adj'ed Be Norm Are p'y Atom %									
MgO	2 0.1s 530.062 366.0 2.61 53.830									
u for tool	3 C 1s284 556 15400 22 728									
y for test	4 7r 3d									
	5 Open the embedded									
	⁶ spreadsheet									
Wel Formula One We	rkback Decimper									
File Edit View Insert E	Elle Edit View Tesert Earent Teser									
	an en 🔊 🗛									
Z4 7.04										
1 XPS Line Ad	j'ed Be Norn A All tools of spreadsheet									
2 01s	530.062 366 available									
3 C1s	284.556 154.									
4 Zr 3d	181.965 159.395 7.33 23.442									
5										
6										
7										
8										
9										
9 10 Areas / Peaks	Composition									

The "Launch Spreadsheet" icon opens the VCI Formula One designer. This provides access to the embedded spreadsheet. A standard set of design tools are available. The redesigned spreadsheet can not be used to change the spreadsheet properties of the ANALYSIS Hawk application. The redesigned spreadsheet can be printed or pasted to other programs.

	VCI	Formula One V	Vorkbook Designer	r	
4	File Ed	lit View Insert	Format Tools		
) 🖻		Cells Row ▶		1
	C	XPS I	Column 🕨 🕨	Width	
r			Sheet 🕨	AutoFit Selection	▲
	1	XPS Line	Printing		
	2	O 1s	Freeze Panes	Default Width 53.830	
ł	3	C 1s	Default Font	154.538 1.00 22.728	
	4	Zr 3d	Object	159.395 7.33 23.442	
F	5		Bring To Front		
	6		Send to Back		
1	7				
	8				
	9				
	10	Aroos & Deel			
Ľ		Areas / Peak	s A composition /		•
Α	djust co	olumn width to fit th	ne selected cells		

The following example displays some of the available tools.

Use the Format > Column > AutoFit Selection to expose all the columns of the spreadsheet.

💵 vci	🐺 VCI Formula One Workbook Designer											
File Edit View Insert Format Tools												
D 🚔 🖶 🚭 🗛 🐰 🛍 🖻 🝼 💆												
C1 XPS Line												
1	XPS	Line	Peak ID	BE Shift	Chem Shift	ådj'ed Be	Norn Area	% Gauss	Sen'y	Group %	Atom %	
2	0 1s			12.000	-1.938	530.062	366.010		2.61	Grp %	53.830	
3	C1s			12.000	044	284.556	154.538		1.00	Grp %	22.728	
4	Zr 3d			12.000	1.965	181.965	159.395		7.33	Grp %	23.442	
5												i l
6												
7												
8												
9												
10	Area	s & Pe	aks A Com	nosition /								•
Ready	7 1000	5 / 1 6		Joomon y								-



The Insert > Column command can be used to enter an empty column.

-														
₩ V	🖤 VCI Formula One Workbook Designer 📃 🔲 🗙													
File	File Edit View Insert Format Tools													
Dı	D 🚔 🖬 🚭 🕼 👗 🖻 🛍 🚿 👧													
	D1 Atom %													
1	2	KPS Line	Atom %	Peak	ID	BE Shift	Chen Shif	t Adj'ed Be	Norm Area	% Gauss	Sen'y	Group %	Atom %	
2	0	D1s	53.829524			12.000	-1.93	8 530.062	366.010		2.61	Grp %	53.830	
3	0	C1s	22.728092			12.000	04	4 284.556	154.538		1.00	Grp %	22.728	
4	Z	(r 3d	23.442384	L		12.000	1.98	5 181.965	159.395		7.33	Grp %	23.442	
5														
	Past	e Special												
	Pasi	te		Pas	ste	special	with							
	6) All		"	.1	aa'' aha	alzad to		Use co	py to	pick u	p		
	C Formulas Values CHECKED IO													
	Atom % column													
Per	ç	Formats				1								
Red		OK	Car	ncel	H	000	40	0	200	_				
						BINDING EN	IERGY - EV							

The Edit > Copy command can be used to copy a complete column and then Edit > Paste Special to fill the empty column. Use the Values mode to insure the numbers are copied and not the formulas.

The rearranged table may be printed or copied to the clipboard. The format tool provides hide column or hide row functions.

The VCI Formula One Designer has been opened with cell protection turned off. This allows the editing of both the cell values and formulas.

DP Viewer

DP Viewer provides an alternate depth profile tool. This tool is in development. Any suggestions for improvement are welcome.

In the first window you may browse to an alternate Scofield table. It is recommended that you use the Programdb.mdb that is in the C:\Program Files\Esca 2000 folder. The Scofield table contained in the Programdb.mdb can be loaded by selecting [Load Table] when the path is showing in the path window.

After the Scofield is loaded the DP Viewer will open. No default depth profile is loaded. Go to the Folder Icon to open a database. You may navigate your network, hard drive or removable media to locate the folder that contains the database of interest. When the database is selected the project/experiment tree for the database will be displayed in the right pane. Find the Raw data or previous Analysis for the depth profile of interest. Select open to load the depth profile data into the DP Viewer.

Locate Scofield	Table		
The Be sure program	path is to ESC ndb.mdb	CA 07	r should be in a Microsoft Jet database file (.mdb). The database should have a
Atomic	floating poi	atomic number, e.g. 79	
TipicalEv	floating poi	binding energy, e.g. 83.98	
CrossSec	floating po	9.58	
Lbe	floating po	low end of range, e.g	
Ube	floating po	low end of range, e.g C	neck the "don't show"
Symol	String	atomic symbol, e.g. A ch	eck box
Name	String	element name, e.g. gold	
Principle	boolean	true if this is strong live	
Please indicate whe	re to find the fi	n the 'Cancel' but cept de	ault values.
C:\Data Folders\ESCA 0)7\programdb.mdb	Browse	
File in edi	it coput in the loaded, unle Use factory provi w this dialog unless error occ	es an error occurs.	Click Load Table



DP Viewer with data loaded. There are no default baselines assigned. Select the spectral region to view. Cycle 1 will be displayed by default.



Select a measurement mode. The Zoom tool is turned off when you enter a measurement mode. The Area mode is shown in this example.

Set the end point average. In the example, 3 points to each side of the selected point will be averaged, for a total of 7 points.

Choose linear or Shirley base line.

To draw a base line click near the spectrum. The end points will jump to the point you click. There is no need to "sneak up" to the current end point, activate it and then drag. When you like the starting point then drag. Now click [Next >] to measure the area.

If this is the first spectrum, for the selected region, then [Is Fit Template] check box will be checked by default. This means the binding energy of the two baseline endpoints can be used to define the base line for the rest of the cycles. Use the [Propagate Reg] button to propagate the area measurement to all cycles of the selected region.

You can move to a different cycle, draw a baseline and select next to make an area measurement. Then check the [Is Fit Template] check box. Now this cycle will become the Template. The original cycle will be automatically unchecked. The new template can be propagated to all cycles.

A cycle can be selected and its endpoints modified. When next is selected the new endpoints will be used to measure the area for this cycle only. Not all cycles need to have the same endpoints.



After the [Next >] button is clicked the area is measured and presented in the table.



The cycle and be changed using the up/down arrows or the mouse wheel when the cycle control is active.



Select a new region using the [Region Control]



To view all cycles at one time select the [Group Region Tab]



Double click the 3D display to bring up the Display Control. May options are available.



The Atom percent is computed and displayed based on the regions/cycles measured. It is only computed for the cycle displayed.



Preview of Depth Profile. This preview can be printed. The chart editor provides font control, background color control and much more.



Exporting Depth Profile to EXCEL

To export the profile to EXCEL select the "Make a Report" button in the DP Viewer toolbar. The "Chart Options Dialog will open. The X axis scale type can be selected as Depth, Etch Time or Cycle. If Depth is selected enter the Etch rate. You may also enter a title for the Depth Chart. Select Continue.

Chart Options D)ialog		
Etch Rate (nm/min)	34.5		
Chart Title	Si02 thickness	s test	
Chart Horiz Axis	Depth	C Etch Time	🔿 Cycle
Edit Chart For	mat Code	Edit Series Format Code	Continue

The Scofield Table locator will be presented if the "don't show" check box is unchecked. Select the [Load Table] button.

🕅 Locate Scofiel	d Table	🔟 Locate Scofield Table							
This program is atte for should be in a Mi have a table named	mpting to locate a scofie icrosoft Jet database file I 'Scofield' with the follow	eld table. The table it is looking (mdb). The database should ring fields:							
Field Name	Туре	Use							
Element	String	Au 4f7							
Atomic	floating point	atomic number, e.g. 79							
TipicalEv	floating point	binding energy, e.g. 83.98							
CrossSec floating point 9.58									
Lbe	low end of range, e.g. 82								
Ube	Jbe floating point low end o								
Symol	String	atomic symbol, e.g. Au							
Name	String	element name, e.g. gold							
Principle	boolean	true if this is strong line							
Please indicate whe default values.	ere to find the file, or push	the 'Cancel' button to accept							
C:\Program Files\Esca	2000\Programdb.mdb	Browse							
File in edit control will be loaded, unless an error occurs.									
	Use factory provided default values Cancel								
☐ Next time, don't show this dialog unless error occurs.									

After the export, the EXCEL application may be minimized but the Indicator will be present in the Windows Status bar at the bottom of the display. If you view the DP Viewer window you will find the following Report Log. It may contain the following Error in the 5 line. This is a result of the security setting in EXCEL for using VBA programming. This is not a fatal error. The Export will be complete. Some special programming features that are not currently activated will not be available. The Error can be eliminated by changing the security setting as described in the Software Installation section of this Instruction.



Trust Center	? 🗙 🤇
Trusted Publishers Trusted Locations Add-ins ActiveX Settings Macro Settings	Macro Settings For macros in documents not in a trusted location: Oisable all macros without notification Isable all macros with notification Oisable all macros except digitally signed macros Enable all macros (not recommended; potentially dangerous code can run)
Message Bar	Developer Macro Settings
External Content Privacy Options	✓ Trust access to the <u>VBA</u> project object model
	OK Cancel

Security setting for EXCEL 2007 used to avoid programming error.



Report with Security setting changed to "Trust" VBA programming.

After the Excel application is loaded with all the depth profile data the Excel Tab will stop blinking in the Windows Application Tray. Select the tab to display the tables.



The first tab will display the depth profile. Atomic percent data tables are summarized in the Summary table. The Composition Table tab displays full data for all cycles.



0	9-	(24 -)∓							M	icrosoft Exc	el
	Home	Insert	Page Layo	ut Form	ulas Da	ta Revi	ew View	/ Develo	oper Ac	robat	
-	🖁 🔏 Cut		Calibri	× 11	· A ·	=		📑 Wrap Te	ext	General	*
Past	Copy		ID 7 II					- ·	Contor a	¢ - 0/	0 .00
*	🔰 🍼 Forma	at Painter	DIU			= = =		merge a	x center *	⊅ ° 70	7 .00 →.0
	Clipboard D Font D Alignment D Number D										
	A1 • (Summary Table for BindingEnergy										
1	Book1										
	А	В	С	D	E	F	G	Н	- I	J	К
1	Summary	Table for B	BindingEne	rgy							
2		Cycle	EtchTime	Depth	Si 2p	O 1s	C 1s				
3	Cycle_1	1	0	0	102.481	532.3162	284.3117				
4	Cycle_2	2	90	51.75	102.8318	532.3728	284.5275				
5	Cycle_3	3	180	103.5	102.7563	532.4178	285.3359				
6	Cycle_4	4	270	155.25	102.9108	532.4284	284.1075				
7	Cycle_5	5	360	207	102.6595	532.4684	286.0585				
8	Cycle_6	6	450	258.75	102.4377	532.4183	283.6938				
9	Cycle_7	7	540	310.5	102.6671	532.388	284.3129				
10	Cycle 8	8	630	362.25	102.4574	532.4308	284.1707				
11	Cycle 9	9	720	414	102.8823	532.3896	284.25				
12	Cycle 10	10	810	465.75	102.3564	532.3679	284.5404				
13	Cycle 11	11	900	517.5	102.5206	532.1038	285.1686				
14	Cycle 12	12	990	569.25	102.5286	532.0783	282.8661				
15	Cycle 13	13	1080	621	102.7439	532.3097	285.8067				
16	Cycle 14	14	1170	672.75	102.4726	532.3386	284.0458				
17	Cycle 15	15	1260	724.5	102.6425	532.4025	285.5477				
18	Cycle 16	16	1350	776.25	102.6945	532.3604	284.1024				
19	Cycle 17	17	1440	828	102.7373	532.3489	283.9015				
20	Cycle 18	18	1530	879.75	103.0859	532.3201	285.9295				
21	Cycle 19	19	1620	931.5	102.7438	532.0421	285.7705				
22	Cycle 20	20	1710	983.25	102.701	532.1041	285.8831				
23	Cycle 21	21	1800	1035	101.9915	531.9995	284.2719				
24	Cycle 22	22	1890	1086.75	98.35308	530.7889	283.8077				
25	Cycle 23	23	1980	1138.5	98.37248	530.8023	284.771				
26	Cycle 24	24	2070	1190.25	98.36657	531.5303	286.417				
27	Cycle 25	25	2160	1242	98.35715	530.6172	286.2361				
28	Cycle 26	26	2250	1293.75	98,36073	531.0209	284.2061				
29	Cycle 27	27	2340	1345.5	98,34097	531,3896	284.9839				
30	Cycle 28	28	2430	1397.25	98,31579	531.049	283.5266				
31	Cycle 29	29	2520	1449	101,5074	531,5988	284.984				
32	_,		2020	2.75							
33	Summarv	Table for A	Area								
34		Cycle	EtchTime	Depth	Si 2p	O 1s	C 1s				
35	Cycle 1	1	0	0	12308.37	52124.95	5397.748				
36	Cycle 2	2	90	51.75	16763.95	67053.01	340.0399				
					20.00000	7.000.01	ition Table	Charte	Charte	(\$ 7	
	At		sa 🖉 Bindin	g Energy	Summary	Compos	scion Tables	Sneet2	sneet3		
Ready	Ready 🔠										

Summary tables in EXCEL

C.	9-	(°I -) -							Mi	crosoft Exc	el	
0	Home	Insert	Page Layo	ut Form	ulas Da	ta Revie	ew View	Develo	per Acr	obat		
	🔏 Cut		Calibri	* 11	A A	= = =	\$\ }~~	📑 Wrap Te	ext	General	-	
Past	e Copy		BIU	-	» - A -			•a• Merge 8	k Center 👻	\$ - %	• €.0 .00	
*	Clipboard	at Painter		Font	5		Alianm	ent	G	Numb	er G	
	A1 • fx Composition table for cycle 1 Decrease Indept (Ctrl=Alt=Shift=Tab)											
Decrease the margin between the												
	BOOKI	D	6	D	F	E.	border	border and the text in the cell.				
1	A	в on table f	or cycle 1	U	E	F	G	н		J	ĸ	
2	Name	Line	Be	Fwhm	Height	Area	CrossSecti	RelArea	Atom%	ScanTime	AdiCs	
3	Si 2p	Si 2p	102.481	2.328425	838	12308.37	0.817	90.99108	31.99923	150	0.9018	
4	0 1s	0 1s	532.3162	2.152805	3666	52124.95	2.93	139.3945	49.02147	150	2.492922	
5	C 1s	C 1s	284.3117	2.162288	341	5397.748	1	53.96842	18.97931	100	1.000168	
6												
7												
8	Composit	ion table f	or cycle 2									
9	Name	Line	Be	Fwhm	Height	Area	CrossSect	RelArea	Atom%	ScanTime	AdjCs	
10	Si 2p	Si 2p	102.8318	2.32059	1078	16763.95	0.817	123.9515	40.41789	150	0.90164	
11	O 1s	O 1s	532.3728	2.102947	4901	67053.01	2.93	179.3231	58.47336	150	2.492819	
12	C 1s	C 1s	284.5275	8.574507	105	340.0399	1	3.400256	1.108749	100	1.000042	
13												
14	Compositi	ion table f	or cuclo 2									
15	Namo	line	or cycle 3	Ewbm	Hoight	Aroa	CrossSoct	PolAroa	Atom%	ScanTimo	AdiCo	
17	Si 2n	si 2n	102 7562	2 254092	1126	16221 20	0.917	110 0252	20 5/100	scannine 150	AUJUS 0.901675	
18	0.1s	0.1s	532 4178	2.234052	4874	67598.06	2 93	119.9352	59 60303	150	2 492737	
19	C 15	C 15	285.3359	8.384813	118	259.491	2.55	2.596023	0.855875	100	0.999571	
20	0.15		20010000	0.001010		2007102		21050020	0.000070	200	0.000000	
21												
22	Compositi	ion table f	or cycle 4									
23	Name	Line	Be	Fwhm	Height	Area	CrossSect	RelArea	Atom%	ScanTime	AdjCs	
24	Si 2p	Si 2p	102.9108	2.420061	1049	16411.27	0.817	121.3487	39.85391	150	0.901604	
25	O 1s	O 1s	532.4284	2.085839	4935	67096.4	2.93	179.4465	58.93466	150	2.492717	
26	C 1s	C 1s	284.1075	9.000133	105	368.9641	1	3.688583	1.211422	100	1.000287	
27												
28												
29	Composit	ion table f	or cycle 5	E h.	11-1-1-1	•	.	Dela	A	0	A direct	
30	Name	Line	Be	Fwhm	Height	Area	CrossSecti	RelArea	Atom%	ScanTime	AdjCs	
31	51 2p	51 2p	102.0595	2.23865	1089	15927.64	0.81/	101 2474	39.025	150	0.901/19	
32	015	015	296 0505	2.083498	4919	07/0/.8	2.93	2 7//170	0.00558	100	0.00015	
33	C 15	0.12	200.0383	7.040398	38	274.1048	1	2.744179	0.505424	100	0.53515	
35												
36	Compositi	ion table f	or cycle 6									
	sempositi	an coord i										

Cycle by Cycle Composition Tables



Depth vs Area Chart.

The data exported to EXCEL allows the Analyst flexibility in organizing, formatting and presenting the data.

Other DP Viewer Functions

Auto-lookup of Scofield crossection. If the label given each region uses the standard XPS transition line format, such as Au 4f, then the Auto-lookup feature will find the Scofield crossection and compute the normalized area.



If the label is not found in the Scofield table then a dialog will be displayed. This dialog may result from errors in entering the label.



In the Capture Application the XPS Line designations can be used to label the Region Names. This provides the correct information to the Auto-Lookup feature in the DP Viewer Application.

Correct entry of XPS Line designations in Depth Profile Dialog will facilitate the auto assignment of Scofield Crossection information.

1 Function Spot Siz Depth Profile ▼ 400 Au on SS ▼ ✓	e Etch Flood Gun Cycles Rotate 00:00:15 off 6 off END	Clear ↑ ↓ × A
Depth Profile Dialog Template Name Au on SS	Create Date Last Modified	× ⊻ ₹
Spot Size Etch 400 00:00:15 Rotate	Flood Gun Cycles Off 6 Set As Default Yes	
RegionName CBE Scar	Region Definitions nned? WW Scans/Time Res# EV/Ster _ a	B C D E
Au 4f 84 1 C 1s 04 1 0 1s 530 530 Fe 710 1	0 00:00:15 4 0 00:00:15 4 0 00:00:15 4 0 00:00:15 4	
Comment	Use XPS line designators a Region Names to enable th Auto-Lookup feature in DI Viewer	as ne P

Line Color Control and Export all Cycles of one Region



Configuring EXCEL Security

Excel allows macro programs to be used to automate the behavior of the Excel application. DP Viewer uses this feature to automate the setup of charts and tables in Excel. There has been increasing levels of security added to Excel releases. The security help to keep dangerious programs from being embedded on computers tied to the internet. Excel XP 2003 and Excel XP 2007 have this protection built in. The protection can be modified so DP Viewer will still operate. The following describes how to change the protection settings.

EXCEL XP 2002

Microsoft Excel - Book1		
Elle Edit Yiew Insert Format □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □	Tools Data Window Help	🍓 Σ + ੈ↓ 🕌 🛍 🦓 100% - 📿 .
A B C 1 2 3 4 5 6	Share Workbook Protection Opline Collaboration Formula Auditing Tools on the Web	G H I
7 8 9 10	Macro ► Customize Options ¥	Macros Alt+F8 Record New Macro Security
12 13 14		Visual Basic Editor Alt+F11 Microsoft Script Editor Alt+Shift+F11

Open EXCEL and make following settings.

OPEN THE SECURITY DIALOG

	Microsoft Exe	el - Book1						Security ?X
	Eile Edit	View Inse	rt F <u>o</u> rmat ∂. ♥ / 牀	Tools Da	ata <u>W</u> indow	Help	$\Sigma \rightarrow \frac{1}{2}$	Security Level Trusted Sources
-	A1	▼ j	£			1 -	2.	C High. Only signed macros from trusted sources will be allowed to run. Unsigned macros are automatically disabled.
1 2 3 4 5 6 7 8 9 10 11 12 13			C	D	E	F	G	Medium. You can choose whether or not to run potentially unsafe macros. Low (not recommended). You are not protected from potentially unsafe macros. Use this setting only if you have virus scanning software installed, or you are sure all documents you open are safe.
14 15 16								OK Cancel

SET SECURITY LEVEL

M	licrosoft E	ксel - Book	:1					Security ?X	
	<u>File E</u> dit	⊻iew I	nsert F <u>o</u> rmat	<u>T</u> ools <u>D</u> ata	<u>W</u> indov	v <u>H</u> elp	- 41	Security Level Trusted Sources	
		5 🖸 (6	\$ _9. ♥ ₺	🖷 🖪 • 🛇	1 2 4	CH + 👹	∑·• Ž↓		ŀ
	AI	•	Ţx.						Ŀ
	Α	В	C	D	E	F	G		L
1									L
2									L
3									L
4									Г
5									E
6									E
7									E
8									E
9								Remove	E
10								Trust all installed add-ins and templates	Г
11								Truck gammistance ded into and templates	Г
12								IV Trust access to visual basic Pr0)ect	E
13								No virus scapper installed	E
14									E
15								OK Cancel	E
16									

SET TRUST SETTING TO ALLOW ACCESS TO VISUAL BASIC PROJECTS.

EXCEL XP 2007

Excel 2007 must be installed with the developer options. With the Developer options selected select the Macro Security dialog. Make settings as shown.

	- (* -) +							Microso	oft Excel
Home	Insert Page La	ayout Fo	ormulas	Data	Review	v View	Developer	Acrobat	
Visual Macros Basic	Record Macro	nces	t Design Mode	Prop View	erties Code Dialog	Source Refresh Data			
A1	Macro Security Customize the macro settings.	security		10.5472			. Anti-		
A 1 2	Press F1 for more	help.	E		F	G	Н	1	1

Trust Center	? 🛛
Trusted Publishers Trusted Locations Add-ins ActiveX Settings Macro Settings Message Bar	Macro Settings For macros in documents not in a trusted location: O Disable all macros without notification Image: Disable all macros with notification Image: Disable all macros with notification Image: Disable all macros except digitally signed macros Image: Disable all macros (not recommended; potentially dangerous code can run) Developer Macro Settings
External Content Privacy Options	✓ Trust access to the <u>V</u> BA project object model OK Cancel

If macros are disabled specify notification is required. Select to trust access to VBA project object model.