

From: Rob Roberts <rob.roberts@bio-logic.us>  
Subject: RE: SP-300 used as Bipotentiostat  
Date: January 31, 2017 at 4:49:47 AM PST  
To: "Brunschiwig, Bruce S." <bsb@caltech.edu>, Bill Eggers  
<bill.eggers@bio-logic.us>, Arnold Forman <arnold.forman@bio-logic.us>  
Cc: "Nunez, Paul D." <pnunez@caltech.edu>, "Fermin, David J."  
<davidcal@caltech.edu>

Bruce,

To run the SP-300 as a bi-stat, it needs to run in CE-to-ground mode which is a selection in the Advanced Settings, and special lead connections. See attached photos of SP-300 connected to RRDE system (two WE, with a common RE and CE). In this mode, the ground leads become the P2 (counter) leads, and the regular P2 (dark blue) leads are not used.

Best regards,

Rob Roberts  
Strategic Marketing and Sales Manager  
Bio-Logic USA  
9050 Executive Park Drive  
Suite 110C  
Knoxville, TN 37923

Tel: 865-769-3800  
Email: rob.roberts@bio-logic.us  
[www.ec-lab.com](http://www.ec-lab.com)

